

AHN myLab® TR-02

Tube Rotator 30 rpm

Usage - ideal for processing of blood samples

Capacity - 16 tubes (4 ESR Tubes and 12 blood collection) tubes with the maximum operational weight capacity of 1.5 kgs

Fixed speed - 30 rpm with an accuracy of +/-2 rpm

Optimal mixing - thanks to gentle circular rotations on the fixed plate at the angle of 38 degrees

Robust and space-saving design - robust metal construction with rubber feet for stable operation / compact design requires minimum of bench space

Easy maintenance and adjustment - the disc plate can be easily removed for cleaning and replacing



Ordering information

Description

AHN myLab® TR-02 Tube Rotator 30 rpm

Cat. No.

7-210-00-0

Delivery package

1 pc. Tube Rotator 30 rpm
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card

Specifications

Run time	Continuous / infinite
Noise level	Noiseless operation
Max. loading capacity	1.5 kg
Protection class	IP 21
Dimensions (WxDxH)	278 x 183 x 252 mm
Weight	1.7 kg

