

Pathogen-proof outer ply

The non-woven spun-bond PP external ply of Schild Line® Professional face masks effectively blocks dust particles and droplets that carry pathogens.

Melt-blown interior layer

Schild Line® Professional face masks feature a melt-blown interior layer that provides superior microfiltration and adsorption of minute microbial pathogens

Non-irritating interior ply

A soft interior ply prevents Schild Line® Professional face masks from irritating the skin while providing a comfortable and near-natural breathing experience

RoHS approved nose pin

All Schild Line® Professional face masks feature an RoHS approved insulated nose pin for tailoring the fit of the mask on the nose bridge to prevent air leakage

Ultrasonic-welded ear loops

The ear loops on Schild Line® Professional face masks reduce pressure on the external ear, providing a comfortable experience thanks to the unique latex-free ultrasonic-welding construction

PROPERTIES

- Top-notch 3-layer face protection masks manufactured from high-quality medicalgrade spun-bond non-woven PP fabric
- Melt-blown interior layer for superior pathogen filtration and user protection
- Soft and gentle interior ply protects the skin from irritation and discomfort
- Comfort enhanced ear loops created with ultrasonic welding technology



TAKE NOTE

- Always wash/sanitize hands before handling your mask
- Change masks as frequently as possible, never exceeding 6 hours with one mask
- Do NOT touch your mask or your face in public places
- Only remove the mask by pulling off the ear loop

FEATURES

- Contaminant-free cleanroom production
- UV sterilized for end-user protection
- Exceptional bacterial filtration efficiency
- Ear loops constructed with innovative ultrasonicwelding technology
- Insulated single core nose pin for tailored fit (RoHS approved)
- Soft and gentle 3-layered mask design
- High-quality, odorless non-woven spun-bond fabric construction
- Free of common irritants (e.g. latex and glass fiber)

ORDERING INFORMATION

Cat. No. 6-500-50-20

Description Schild Line® Professiona

Face Masks

Pack 50 x 2 x 25 pcs. = 2500 pcs.

