

ahn myBottle

Carboys

Range of Applications - from dispensing DI water to mixing reagents, bulk vaccines and media

Secure - leakproof design ensures safe transport

Round bottom - ensures even mixing and curved inner corners for easy cleaning

Premium Materials - made from USFDA approved raw material for high quality, durability and safety / available in PP, HDPE, and LDPE, ranging from 10L to 50L, with both sterile and non-sterile options

Versatile - ideal for storage, mixing, dispensing, and transporting liquids

User-Friendly - ergonomic handles, wide-mouth openings, and spigot for easy use

Options - various materials and capacities

Precise Identification - unique lot numbers and laser printed graduations


Certified USP Class VI material and animal component free

Compliant to ISO 10993/TSE/BSE Free and BPA Free

Profiles - Extractable and Leachable (E&L) Profile available

Sterility - Available with sterility assurance level 10⁻⁶ as per ISO 11137

Optional with unique serial number - for easy identification based on customer requirement



Printed graduations



Carboy spigot

Ordering information

Carboys, HDPE

Description	Sterile	Cat. No.
Carboy 10L, w/o cap		10-035-01-0
Carboy 15L, w/o cap		10-036-01-0
Carboy 20L, w/o cap		10-037-01-0

Carboys, HDPE

Description	Sterile	Cat. No.
Carboy 25L, w/o cap		10-038-01-0
Carboy 50L, w/o cap		10-039-01-0
Carboy 10L, w/ cap		10-040-01-0

Carboys, HDPE

Description	Sterile	Cat. No.
Carboy 15L, w/ cap		10-041-01-0
Carboy 20L, w/ cap		10-042-01-0
Carboy 25L, w/ cap		10-043-01-0
Carboy 50L, w/ cap		10-044-01-0
Carboy 10L, w/ spigot, w/o cap		10-045-01-0
Carboy 15L, w/ spigot, w/o cap		10-046-01-0
Carboy 20L, w/ spigot, w/o cap		10-047-01-0
Carboy 25L, w/ spigot, w/o cap		10-048-01-0
Carboy 50L, w/ spigot, w/o cap		10-049-01-0
Carboy 10L, w/ spigot, w/ cap		10-050-01-0
Carboy 15L, w/ spigot, w/ cap		10-051-01-0
Carboy 20L, w/ spigot, w/ cap		10-052-01-0
Carboy 25L, w/ spigot, w/ cap		10-053-01-0
Carboy 50L, w/ spigot, w/ cap		10-054-01-0

Carboys, HDPE, sterile

Description	Sterile	Cat. No.
Carboy 10L, w/o cap	■	10-055-01-0
Carboy 15L, w/o cap	■	10-056-01-0
Carboy 20L, w/o cap	■	10-057-01-0
Carboy 25L, w/o cap	■	10-058-01-0
Carboy 50L, w/o cap	■	10-059-01-0
Carboy 10L, w/ cap	■	10-060-01-0
Carboy 15L, w/ cap	■	10-061-01-0
Carboy 20L, w/ cap	■	10-062-01-0
Carboy 25L, w/ cap	■	10-063-01-0
Carboy 50L, w/ cap	■	10-064-01-0
Carboy 10L, w/ spigot, w/o cap	■	10-065-01-0
Carboy 15L, w/ spigot, w/o cap	■	10-066-01-0
Carboy 20L, w/ spigot, w/o cap	■	10-067-01-0
Carboy 25L, w/ spigot, w/o cap	■	10-068-01-0
Carboy 50L, w/ spigot, w/o cap	■	10-069-01-0

Carboys, LDPE

Description	Sterile	Cat. No.
Carboy 10L, w/o cap		10-070-01-0
Carboy 15L, w/o cap		10-071-01-0

Carboys, LDPE

Description	Sterile	Cat. No.
Carboy 20L, w/o cap		10-072-01-0
Carboy 25L, w/o cap		10-073-01-0
Carboy 50L, w/o cap		10-074-01-0
Carboy 10L, w/ cap		10-075-01-0
Carboy 15L, w/ cap		10-076-01-0
Carboy 20L, w/ cap		10-077-01-0
Carboy 25L, w/ cap		10-078-01-0
Carboy 50L, w/ cap		10-079-01-0
Carboy 10L, w/ spigot, w/o cap		10-080-01-0
Carboy 15L, w/ spigot, w/o cap		10-081-01-0
Carboy 20L, w/ spigot, w/o cap		10-082-01-0
Carboy 25L, w/ spigot, w/o cap		10-083-01-0
Carboy 50L, w/ spigot, w/o cap		10-084-01-0
Carboy 10L, w/ spigot, w/ cap		10-085-01-0
Carboy 15L, w/ spigot, w/ cap		10-086-01-0
Carboy 20L, w/ spigot, w/ cap		10-087-01-0
Carboy 25L, w/ spigot, w/ cap		10-088-01-0
Carboy 50L, w/ spigot, w/ cap		10-089-01-0

Carboys, LDPE, sterile

Description	Sterile	Cat. No.
Carboy 10L, w/o cap,	■	10-090-01-0
Carboy 15L, w/o cap	■	10-091-01-0
Carboy 20L, w/o cap	■	10-092-01-0
Carboy 25L, w/o cap	■	10-093-01-0
Carboy 50L, w/o cap	■	10-094-01-0
Carboy 10L, w/ cap	■	10-095-01-0
Carboy 15L, w/ cap	■	10-096-01-0
Carboy 20L, w/ cap	■	10-097-01-0
Carboy 25L, w/ cap	■	10-098-01-0
Carboy 50L, w/ cap	■	10-099-01-0
Carboy 10L, w/ spigot, w/o cap	■	10-100-01-0
Carboy 15L, w/ spigot, w/o cap	■	10-101-01-0
Carboy 20L, w/ spigot, w/o cap	■	10-102-01-0
Carboy 25L, w/ spigot, w/o cap	■	10-103-01-0
Carboy 50L, w/ spigot, w/o cap	■	10-104-01-0



Carboys, PP

Description	Sterile	Cat. No.
Carboy 10L, w/o cap		10-000-01-0
Carboy 15L, w/o cap		10-001-01-0
Carboy 20L, w/o cap		10-002-01-0
Carboy 25L, w/o cap		10-003-01-0
Carboy 50L, w/o cap		10-004-01-0
Carboy 10L, w/ cap		10-005-01-0
Carboy 15L, w/ cap		10-006-01-0
Carboy 20L, w/ cap		10-007-01-0
Carboy 25L, w/ cap		10-008-01-0
Carboy 50L, w/ cap		10-009-01-0
Carboy 10L, w/ spigot, w/o cap		10-010-01-0
Carboy 15L, w/ spigot, w/o cap		10-011-01-0
Carboy 20L, w/ spigot, w/o cap		10-012-01-0
Carboy 25L, w/ spigot, w/o cap		10-013-01-0
Carboy 50L, w/ spigot, w/o cap		10-014-01-0
Carboy 10L, w/ spigot, w/ cap		10-015-01-0
Carboy 15L, w/ spigot, w/ cap		10-016-01-0
Carboy 20L, w/ spigot, w/ cap		10-017-01-0
Carboy 25L, w/ spigot, w/ cap		10-018-01-0
Carboy 50L, w/ spigot, w/ cap		10-019-01-0

Carboys, PP, sterile

Description	Sterile	Cat. No.
Carboy 10L, w/o cap	■	10-020-01-0
Carboy 15L, w/o cap	■	10-021-01-0
Carboy 20L, w/o cap	■	10-022-01-0
Carboy 25L, w/o cap	■	10-023-01-0
Carboy 50L, w/o cap	■	10-024-01-0
Carboy 10L, w/ cap	■	10-025-01-0
Carboy 15L, w/ cap	■	10-026-01-0
Carboy 20L, w/ cap	■	10-027-01-0
Carboy 25L, w/ cap	■	10-028-01-0
Carboy 50L, w/ cap	■	10-029-01-0
Carboy 10L, w/ spigot, w/o cap	■	10-030-01-0
Carboy 15L, w/ spigot, w/o cap	■	10-031-01-0
Carboy 20L, w/ spigot, w/o cap	■	10-032-01-0
Carboy 25L, w/ spigot, w/o cap	■	10-033-01-0
Carboy 50L, w/ spigot, w/o cap	■	10-034-01-0



myBottle Screw Caps

Description	Pack Type	Sterile	Cat. No.
Screw cap for Carboy 83B PP			10-900-01-0
Screw cap for Carboy 83B PP	bulk		10-900-05-0
Screw cap for Carboy 83B PP		■	10-901-01-0
Screw cap for Carboy 83B PP	bulk	■	10-901-05-0

