

ahn myBottle

Media Bottles



Volume – square shaped, in 30 mL to 2 litres, with graduated markings

Design - high-quality PETG/PET materials, minimising additives and leachables. Suits biopharmaceutical processing, critical reagent storage

Security - shatterproof/leakproof for sample integrity / tamper-proof seals

Easy-to-clean - smooth surface prevents residue build-up

Manufactured in Class 8 cleanroom from USP Class VI Raw material

Temperature resistant - from 60°C to -70°C

Filter/plasticiser-free - uncompromised sample integrity

Versatile - Gamma sterile, non-pyrogenic, non-cytotoxic. TSE/BSE Compliant and BPA Free material. Compliant to ISO 10993

PC Media Bottles are available on request

Ordering information

Media Bottles, PETG

Volume	Sterile	Pack Type	Sales Unit	Cat. No.
Media Bottle w/ cap, 30 mL		bag	case / 12x24	AHN-PETG-30-B
Media Bottle w/ cap, 30 mL	■	bag	case / 4x24	AHN-PETG-30-S
Media Bottle w/ cap, 30 mL	■	bag	case / 12x24	AHN-PETG-30-S-B
Media Bottle w/ cap, 60 mL		bag	case / 8x24	AHN-PETG-60-B
Media Bottle w/ cap, 60 mL	■	bag	case / 4x24	AHN-PETG-60-S
Media Bottle w/ cap, 60 mL	■	bag	case / 8x24	AHN-PETG-60-S-B
Media Bottle w/ cap, 125 mL		bag	case / 4x24	AHN-PETG-125-B
Media Bottle w/ cap, 125 mL	■	bag	case / 4x24	AHN-PETG-125-S
Media Bottle w/ cap, 125 mL	■	bag	case / 4x24	AHN-PETG-125-S-B
Media Bottle w/ cap, 250 mL		bag	case / 2x30	AHN-PETG-250-B
Media Bottle w/ cap, 250 mL	■	bag	case / 2x24	AHN-PETG-250-S
Media Bottle w/ cap, 250 mL	■	bag	case / 2x30	AHN-PETG-250-S-B
Media Bottle w/ cap, 500 mL		bag	case / 2x24	AHN-PETG-500-B
Media Bottle w/ cap, 500 mL	■	bag	case / 2x24	AHN-PETG-500-S
Media Bottle w/ cap, 500 mL	■	bag	case / 2x24	AHN-PETG-500-S-B
Media Bottle w/ cap, 1L		bag	case / 2x12	AHN-PETG-1000-B
Media Bottle w/ cap, 1L	■	bag	case / 2x12	AHN-PETG-1000-S
Media Bottle w/ cap, 2L		bag	case / 2x6	AHN-PETG-2000-B
Media Bottle w/ cap, 2L	■	bag	case / 2x6	AHN-PETG-2000-S

Media Bottles, PET

Volume	Sterile	Pack Type	Sales Unit	Cat. No.
Media Bottle w/ cap, 125 mL		bag	case / 4x24	AHN-PET-125
Media Bottle w/ cap, 125 mL	■	bag	case / 4x24	AHN-PET-125-S
Media Bottle w/ cap, 250 mL		bag	case / 2x30	AHN-PET-250
Media Bottle w/ cap, 250 mL	■	bag	case / 2x30	AHN-PET-250-S
Media Bottle w/ cap, 500 mL		bag	case / 2x24	AHN-PET-500
Media Bottle w/ cap, 500 mL	■	bag	case / 2x24	AHN-PET-500-S
Media Bottle w/ cap, 1L		bag	case / 2x12	AHN-PET-1000
Media Bottle w/ cap, 1L	■	bag	case / 2x12	AHN-PET-1000-S