

ahn myBottle

Reagent Bottles

Material - US FDA approved raw material of USP class VI, low leachable and extractable. High density polyethylene (HDPE), low density polyethylene (LDPE) and polypropylene (PP) configurations. for chemical resistance and thermal stability

Caps - All caps are made of polypropylene material

Colour - translucent and amber colour variants

Freezer-friendly LDPE / HDPE can withstand temperatures from -100 to 120°C.

PP versions are heat-stable 0 to 135°C, making them perfect for autoclaving

Sizes - Narrow mouth Bottles sizes include 4 mL, 8 mL, 15 mL, 30 mL, 60 mL, 125 mL, 250 mL, 500 mL, 1L, 2L and 4L sizes. Wide mouth bottles come in 30 mL, 60 mL, 125 mL, 175 mL, 250 mL, 500 mL, 1L, 2L and 4L sizes

Chemically resistant/leakproof - suited to both short and long-term storage

Manufacturing - produced in a Class 100K cleanroom facility



Ordering information

Reagent Bottles, HDPE

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 4 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-4
Bottle w/ cap, 4 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-4-A
Bottle w/ cap, 4 mL	narrow	round	translucent		bag	case / 2000	AHN-HD-N-4-B
Bottle w/ cap, 4 mL	narrow	round	amber		bag	case / 2000	AHN-HD-N-4-B-A
Bottle w/ cap, 4 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-4-S
Bottle w/ cap, 4 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-4-S-A
Bottle w/ cap, 8 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-8
Bottle w/ cap, 8 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-8-A
Bottle w/ cap, 8 mL	narrow	round	translucent		bag	case / 2000	AHN-HD-N-8-B
Bottle w/ cap, 8 mL	narrow	round	amber		bag	case / 2000	AHN-HD-N-8-B-A
Bottle w/ cap, 8 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-8-S
Bottle w/ cap, 8 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-8-S-A
Bottle w/ cap, 15 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-15
Bottle w/ cap, 15 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-15-A
Bottle w/ cap, 15 mL	narrow	round	translucent		bag	case / 2000	AHN-HD-N-15-B
Bottle w/ cap, 15 mL	narrow	round	amber		bag	case / 2000	AHN-HD-N-15-B-A
Bottle w/ cap, 15 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-15-S
Bottle w/ cap, 15 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-15-S-A



+49(0)3631/65242-0



info@ahn-bio.com



www.ahn-bio.com



Reagent Bottles, HDPE

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 30 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-30
Bottle w/ cap, 30 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-30-A
Bottle w/ cap, 30 mL	narrow	round	translucent		bag	case / 1000	AHN-HD-N-30-B
Bottle w/ cap, 30 mL	narrow	round	amber		bag	case / 1000	AHN-HD-N-30-B-A
Bottle w/ cap, 30 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-30-S
Bottle w/ cap, 30 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-30-S-A
Bottle w/ cap, 60 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-60
Bottle w/ cap, 60 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-60-A
Bottle w/ cap, 60 mL	narrow	round	translucent		bag	case / 1000	AHN-HD-N-60-B
Bottle w/ cap, 60 mL	narrow	round	amber		bag	case / 1000	AHN-HD-N-60-B-A
Bottle w/ cap, 60 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-60-S
Bottle w/ cap, 60 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-60-S-A
Bottle w/ cap, 125 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-125
Bottle w/ cap, 125 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-125-A
Bottle w/ cap, 125 mL	narrow	round	translucent		bag	case / 500	AHN-HD-N-125-B
Bottle w/ cap, 125 mL	narrow	round	amber		bag	case / 500	AHN-HD-N-125-B-A
Bottle w/ cap, 125 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-125-S
Bottle w/ cap, 125 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-125-S-A
Bottle w/ cap, 250 mL	narrow	round	translucent		bag	case / 6x12	AHN-HD-N-250
Bottle w/ cap, 250 mL	narrow	round	amber		bag	case / 6x12	AHN-HD-N-250-A
Bottle w/ cap, 250 mL	narrow	round	translucent		bag	case / 250	AHN-HD-N-250-B
Bottle w/ cap, 250 mL	narrow	round	amber		bag	case / 250	AHN-HD-N-250-B-A
Bottle w/ cap, 250 mL	narrow	round	translucent	■	bag	case / 6x12	AHN-HD-N-250-S
Bottle w/ cap, 250 mL	narrow	round	amber	■	bag	case / 6x12	AHN-HD-N-250-S-A
Bottle w/ cap, 500 mL	narrow	round	translucent		bag	case / 4x12	AHN-HD-N-500
Bottle w/ cap, 500 mL	narrow	round	amber		bag	case / 4x12	AHN-HD-N-500-A
Bottle w/ cap, 500 mL	narrow	round	translucent		bag	case / 125	AHN-HD-N-500-B
Bottle w/ cap, 500 mL	narrow	round	amber		bag	case / 125	AHN-HD-N-500-B-A
Bottle w/ cap, 500 mL	narrow	round	translucent	■	bag	case / 4x12	AHN-HD-N-500-S
Bottle w/ cap, 500 mL	narrow	round	amber	■	bag	case / 4x12	AHN-HD-N-500-S-A
Bottle w/ cap, 1L	narrow	round	translucent		bag	case / 2x12	AHN-HD-N-1000
Bottle w/ cap, 1L	narrow	round	amber		bag	case / 2x12	AHN-HD-N-1000-A
Bottle w/ cap, 1L, sep. caps	narrow	round	translucent		bag	case / 50	AHN-HD-N-1000-SC
Bottle w/ cap, 1L, sep. caps	narrow	round	amber		bag	case / 50	AHN-HD-N-1000-SC-A
Bottle w/ cap, 2L	narrow	round	translucent		bag	case / 1x6	AHN-HD-N-2000
Bottle w/ cap, 4L	narrow	round	translucent		bag	case / 1x6	AHN-HD-N-4000
Bottle w/ cap, 30 mL	wide	round	translucent		bag	case / 6x12	AHN-HD-W-30
Bottle w/ cap, 30 mL	wide	round	amber		bag	case / 6x12	AHN-HD-W-30-A
Bottle w/ cap, 30 mL	wide	round	translucent		bag	case / 1000	AHN-HD-W-30-B
Bottle w/ cap, 30 mL	wide	round	amber		bag	case / 1000	AHN-HD-W-30-B-A
Bottle w/ cap, 30 mL	wide	round	translucent	■	bag	case / 6x12	AHN-HD-W-30-S

Reagent Bottles, HDPE

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 30 mL	wide	round	amber	■	bag	case / 6x12	AHN-HD-W-30-S-A
Bottle w/ cap, 30 mL , individ	wide	round	translucent	■	bag	case / 500	AHN-HD-W-30-S-IV
Bottle w/ cap, 30 mL , individ	wide	round	amber	■	bag	case / 500	AHN-HD-W-30-S-IV-A
Bottle w/ cap, 60 mL	wide	round	translucent		bag	case / 6x12	AHN-HD-W-60
Bottle w/ cap, 60 mL	wide	round	amber		bag	case / 6x12	AHN-HD-W-60-A
Bottle w/ cap, 60 mL	wide	round	translucent		bag	case / 1000	AHN-HD-W-60-B
Bottle w/ cap, 60 mL	wide	round	amber		bag	case / 1000	AHN-HD-W-60-B-A
Bottle w/ cap, 60 mL	wide	round	translucent	■	bag	case / 6x12	AHN-HD-W-60-S
Bottle w/ cap, 60 mL	wide	round	amber	■	bag	case / 6x12	AHN-HD-W-60-S-A
Bottle w/ cap, 60 mL , individ	wide	round	translucent	■	bag	case / 250	AHN-HD-W-60-S-IV
Bottle w/ cap, 60 mL , individ	wide	round	amber	■	bag	case / 250	AHN-HD-W-60-S-IV-A
Bottle w/ cap, 125 mL	wide	round	translucent		bag	case / 6x12	AHN-HD-W-125
Bottle w/ cap, 125 mL	wide	round	amber		bag	case / 6x12	AHN-HD-W-125-A
Bottle w/ cap, 125 mL	wide	round	translucent		bag	case / 500	AHN-HD-W-125-B
Bottle w/ cap, 125 mL	wide	round	amber		bag	case / 500	AHN-HD-W-125-B-A
Bottle w/ cap, 125 mL	wide	round	translucent	■	bag	case / 6x12	AHN-HD-W-125-S
Bottle w/ cap, 125 mL	wide	round	amber	■	bag	case / 6x12	AHN-HD-W-125-S-A
Bottle w/ cap, 125 mL , individ	wide	round	translucent	■	bag	case / 125	AHN-HD-W-125-S-IV
Bottle w/ cap, 125 mL , individ	wide	round	amber	■	bag	case / 125	AHN-HD-W-125-S-IV-A
Bottle w/ cap, 250 mL	wide	round	translucent		bag	case / 6x12	AHN-HD-W-250
Bottle w/ cap, 250 mL	wide	round	amber		bag	case / 6x12	AHN-HD-W-250-A
Bottle w/ cap, 250 mL	wide	round	translucent		bag	case / 250	AHN-HD-W-250-B
Bottle w/ cap, 250 mL	wide	round	amber		bag	case / 250	AHN-HD-W-250-B-A
Bottle w/ cap, 250 mL	wide	round	translucent	■	bag	case / 6x12	AHN-HD-W-250-S
Bottle w/ cap, 250 mL	wide	round	amber	■	bag	case / 6x12	AHN-HD-W-250-S-A
Bottle w/ cap, 250 mL , individ	wide	round	translucent	■	bag	case / 60	AHN-HD-W-250-S-IV
Bottle w/ cap, 250 mL , individ	wide	round	amber	■	bag	case / 60	AHN-HD-W-250-S-IV-A
Bottle w/ cap, 500 mL	wide	round	translucent		bag	case / 4x12	AHN-HD-W-500
Bottle w/ cap, 500 mL	wide	round	amber		bag	case / 4x12	AHN-HD-W-500-A
Bottle w/ cap, 500 mL	wide	round	translucent		bag	case / 125	AHN-HD-W-500-B
Bottle w/ cap, 500 mL	wide	round	amber		bag	case / 125	AHN-HD-W-500-B-A
Bottle w/ cap, 500 mL	wide	round	translucent	■	bag	case / 4x12	AHN-HD-W-500-S
Bottle w/ cap, 500 mL	wide	round	amber	■	bag	case / 4x12	AHN-HD-W-500-S-A
Bottle w/ cap, 500 mL , individ	wide	round	translucent	■	bag	case / 30	AHN-HD-W-500-S-IV
Bottle w/ cap, 500 mL , individ	wide	round	amber	■	bag	case / 30	AHN-HD-W-500-S-IV-A
Bottle w/ cap, 1L	wide	round	translucent		bag	case / 2x12	AHN-HD-W-1000
Bottle w/ cap, 1L	wide	round	amber		bag	case / 2x12	AHN-HD-W-1000-A
Bottle w/ cap, 1L , sep. caps	wide	round	translucent		bag	case / 50	AHN-HD-W-1000-SC
Bottle w/ cap, 1L , sep. caps	wide	round	amber		bag	case / 50	AHN-HD-W-1000-SC-A
Bottle w/ cap, 2L	wide	round	translucent		bag	case / 1x6	AHN-HD-W-2000
Bottle w/ cap, 2L	wide	round	amber		bag	case / 1x6	AHN-HD-W-2000-A



Reagent Bottles, HDPE

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 4L	wide	round	translucent		bag	case / 1x6	AHN-HD-W-4000
Bottle w/ cap, 4L	wide	round	amber		bag	case / 1x6	AHN-HD-W-4000-A
Bottle w/ cap, 60 mL	wide	square	translucent		bag	case / 6 x 12	AHN-HD-SQ-60
Bottle w/ cap, 60 mL	wide	square	amber		bag	case / 6 x 12	AHN-HD-SQ-60-A
Bottle w/ cap, 175 mL	wide	square	translucent		bag	case / 6 x 12	AHN-HD-SQ-175
Bottle w/ cap, 175 mL	wide	square	amber		bag	case / 6 x 12	AHN-HD-SQ-175-A
Bottle w/ cap, 250 mL	wide	square	translucent		bag	case / 6 x 12	AHN-HD-SQ-250
Bottle w/ cap, 250 mL	wide	square	amber		bag	case / 6 x 12	AHN-HD-SQ-250-A
Bottle w/ cap, 500 mL	wide	square	translucent		bag	case / 4 x 12	AHN-HD-SQ-500
Bottle w/ cap, 500 mL	wide	square	amber		bag	case / 4 x 12	AHN-HD-SQ-500-A
Bottle w/ cap, 1L	wide	square	translucent		bag	case / 2 x 12	AHN-HD-SQ-1000
Bottle w/ cap, 1L	wide	square	amber		bag	case / 2 x 12	AHN-HD-SQ-1000-A
Bottle w/ cap, 2L w/ handle	wide	square	translucent		bag	case / 1 x 6	AHN-HD-SQ-2000
Bottle w/ cap, 4L, w/ handle	wide	square	translucent		bag	case / 1 x 6	AHN-HD-SQ-4000
Bottle w/ cap, 4L, w/ handle	wide	round	translucent		bag	case / 1 x 6	AHN-HD-RB-4000

Reagent Bottles, PP

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 4 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-4
Bottle w/ cap, 4 mL	narrow	round	translucent		bag	case / 2000	AHN-PP-N-4-B
Bottle w/ cap, 8 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-8
Bottle w/ cap, 8 mL	narrow	round	translucent		bag	case / 2000	AHN-PP-N-8-B
Bottle w/ cap, 15 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-15
Bottle w/ cap, 15 mL	narrow	round	translucent		bag	case / 2000	AHN-PP-N-15-B
Bottle w/ cap, 30 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-30
Bottle w/ cap, 30 mL	narrow	round	translucent		bag	case / 1000	AHN-PP-N-30-B
Bottle w/ cap, 60 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-60
Bottle w/ cap, 60 mL	narrow	round	translucent		bag	case / 1000	AHN-PP-N-60-B
Bottle w/ cap, 125 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-125
Bottle w/ cap, 125 mL	narrow	round	translucent		bag	case / 500	AHN-PP-N-125-B
Bottle w/ cap, 250 mL	narrow	round	translucent		bag	case / 6x12	AHN-PP-N-250
Bottle w/ cap, 250 mL	narrow	round	translucent		bag	case / 250	AHN-PP-N-250-B
Bottle w/ cap, 500 mL	narrow	round	translucent		bag	case / 4x12	AHN-PP-N-500
Bottle w/ cap, 500 mL	narrow	round	translucent		bag	case / 125	AHN-PP-N-500-B
Bottle w/ cap, 1L	narrow	round	translucent		bag	case / 2x12	AHN-PP-N-1000
Bottle w/ cap, 1L, sep. caps	narrow	round	translucent		bag	case / 50	AHN-PP-N-1000-SC
Bottle w/ cap, 2L	narrow	round	translucent		bag	case / 1 x 6	AHN-PP-N-2000
Bottle w/ cap, 4L	narrow	round	translucent		bag	case / 1 x 6	AHN-PP-N-4000
Bottle w/ cap, 30 mL	wide	round	translucent		bag	case / 6x12	AHN-PP-W-30
Bottle w/ cap, 30 mL	wide	round	translucent		bag	case / 1000	AHN-PP-W-30-B
Bottle w/ cap, 60 mL	wide	round	translucent		bag	case / 6x12	AHN-PP-W-60
Bottle w/ cap, 60 mL	wide	round	translucent		bag	case / 1000	AHN-PP-W-60-B
Bottle w/ cap, 125 mL	wide	round	translucent		bag	case / 6x12	AHN-PP-W-125
Bottle w/ cap, 125 mL	wide	round	translucent		bag	case / 500	AHN-PP-W-125-B

Reagent Bottles, PP

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 250 mL	wide	round	translucent		bag	case / 6x12	AHN-PP-W-250
Bottle w/ cap, 250 mL	wide	round	translucent		bag	case / 250	AHN-PP-W-250-B
Bottle w/ cap, 500 mL	wide	round	translucent		bag	case / 4x12	AHN-PP-W-500
Bottle w/ cap, 500 mL	wide	round	translucent		bag	case / 125	AHN-PP-W-500-B
Bottle w/ cap, 1L	wide	round	translucent		bag	case / 2x12	AHN-PP-W-1000
Bottle w/ cap, 1L , sep. caps	wide	round	translucent		bag	case / 50	AHN-PP-W-1000-SC
Bottle w/ cap, 2L	wide	round	translucent		bag	case / 1x6	AHN-PP-W-2000
Bottle w/ cap, 4L	wide	round	translucent		bag	case / 1x6	AHN-PP-W-4000
Bottle w/ cap, 1L , heavy duty	wide	round	translucent		bag	case / 1x6	AHN-PP-H-1000
Bottle w/ cap, 2L , heavy duty	wide	round	translucent		bag	case / 1x6	AHN-PP-H-2000
Bottle w/ cap, 4L , heavy duty	wide	round	translucent		bag	case / 1x6	AHN-PP-H-4000
Bottle w/ cap, 60 mL	wide	square	translucent		bag	case / 6x12	AHN-PP-SQ-60
Bottle w/ cap, 175 mL	wide	square	translucent		bag	case / 6x12	AHN-PP-SQ-175
Bottle w/ cap, 250 mL	wide	square	translucent		bag	case / 6x12	AHN-PP-SQ-250
Bottle w/ cap, 500 mL	wide	square	translucent		bag	case / 4x12	AHN-PP-SQ-500
Bottle w/ cap, 1L	wide	square	translucent		bag	case / 2x12	AHN-PP-SQ-1000
Bottle w/ cap, 2L w/ handle	wide	square	translucent		bag	case / 1x6	AHN-PP-SQ-2000
Bottle w/ cap, 4L w/ handle	wide	square	translucent		bag	case / 1x6	AHN-PP-SQ-4000
Bottle w/ cap, 4L , w/ handle	wide	round	translucent		bag	case / 1x6	AHN-PP-RB-4000

Reagent Bottles, LPDE

Description	Neck	Bottom	Colour	Sterile	Pack Type	Sales Unit	Cat. No.
Bottle w/ cap, 4 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-4
Bottle w/ cap, 8 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-8
Bottle w/ cap, 15 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-15
Bottle w/ cap, 30 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-30
Bottle w/ cap, 60 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-60
Bottle w/ cap, 125 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-125
Bottle w/ cap, 250 mL	narrow	round	translucent		bag	case / 6x12	AHN-LD-N-250
Bottle w/ cap, 500 mL	narrow	round	translucent		bag	case / 4x12	AHN-LD-N-500
Bottle w/ cap, 1L	narrow	round	translucent		bag	case / 2x12	AHN-LD-N-1000
Bottle w/ cap, 30 mL	wide	round	translucent		bag	case / 6x12	AHN-LD-W-30
Bottle w/ cap, 60 mL	wide	round	translucent		bag	case / 6x12	AHN-LD-W-60
Bottle w/ cap, 125 mL	wide	round	translucent		bag	case / 6x12	AHN-LD-W-125
Bottle w/ cap, 250 mL	wide	round	translucent		bag	case / 6x12	AHN-LD-W-250
Bottle w/ cap, 500 mL	wide	round	translucent		bag	case / 4x12	AHN-LD-W-500
Bottle w/ cap, 1L	wide	round	translucent		bag	case / 2x12	AHN-LD-W-1000
Bottle w/ cap, 2L	narrow	round	translucent		bag	case / 1 x 6	AHN-LD-N-2000
Bottle w/ cap, 4L	narrow	round	translucent		bag	case / 1 x 6	AHN-LD-N-4000
Bottle w/ cap, 2L	wide	round	translucent		bag	case / 1 x 6	AHN-LD-W-2000
Bottle w/ cap, 4L	wide	round	translucent		bag	case / 1 x 6	AHN-LD-W-4000

