

ahn myLab[®]

Single Plate Magnetic Stirrer 1500 rpm SM-01

Usage - in laboratories, pharmaceutical industries, chemical industries

Maximum speed - from 15 to 1500 rpm with an interval of 50 rpm / software-driven gradual acceleration

Time setting - ranging from 1 - 99 minutes / continuous operation feature

Improved mixing - pulse mode where the mixing direction changes every 30 seconds

Maintenance-free design - water spill, corrosion and chemically resistant / includes no moving parts / high-quality stirring bar

User-friendly interface and design - digital display showing the real-time speed and the remaining time / last run memory function / extremely thin (only 16.0 mm thick) / flat surface enables easy wiping and cleaning



Ordering information

Description	Cat. No.
myLab [®] SM-01 Single Plate Magnetic Stirrer 1500 rpm	7-300-00-0

Optional accessories	Cat. No.
Magnetic stirring bar - PTFE coated, 25x8 mm	CRS-MAG

Delivery Package
1 pc. Magnetic Stirrer 1500 rpm
1 pc. Stir bar (PTFE coated, 25x8 mm)
1 pc. Instruction manual
1 pc. Warranty card

Specifications	
Max. stirring quantity	800 mL
Recommended stir bar length	25 mm
Protection class	IP65
Dimensions (WxDxH)	156x110x13.2 mm
Weight	400 g

