

ahn myPette® EL³


EFFORTLESS

PIPETTING
TIP LOADING
TIP EJECTION



Effortless pipetting

optimal weight distribution and lowest plunger force requirement through innovative spring and seal design



Effortless tip ejection

to ensure continuous uninterrupted pipetting without strain tip ejection



Magnet-assisted piston

for consistently precise results




Calibrated

as per ISO 8655-6



UV resistant



Effortless tip loading

for a reliable and consistent seal between the pipette and the tip without applying excessive force



Colour coding
for easy identification



Volume lock
prevents accidental changes in the volume and guarantees constant pipetting



Large 4-digit display
Precise setting of the volume through a 4-digit and easy-to-read display



Highly durable material
to protect against chemical corrosion and physical stresses



Fully autoclavable



Single Channel Pipettes

Fixed Volume

Cat. No.	Volume Range (µL)	Increments (µL)	Colour Code
8-000-32-9	2.5	-	●
8-001-32-9	5	-	●
8-002-32-9	10	-	●
8-003-32-9	20	-	●
8-004-32-9	25	-	●
8-007-32-9	50	-	●
8-008-32-9	100	-	●
8-009-32-9	200	-	●
8-010-32-9	250	-	●
8-011-32-9	500	-	●
8-012-32-9	1000	-	●

Ordering Information

Variable Volume

Cat. No.	Volume Range (µL)	Increments (µL)	Colour Code
8-100-32-9	0.1-2.5	0.002	●
8-101-32-9	0.5-10	0.02	●
8-102-32-9	2-20	0.02	●
8-103-32-9	5-50	0.1	●
8-104-32-9	10-100	0.1	●
8-105-32-9	20-200	0.2	●
8-106-32-9	100-1000	1	●
8-107-32-9	500-5000	10	●
8-108-32-9	1000-10000	20	●


Multi Channel Pipettes

8 Channel

Cat. No.	Volume Range (µL)	Increments (µL)	Colour Code
8-201-32-9	0.5-10	0.02	●
8-205-32-9	2-20	0.02	●
8-202-32-9	5-50	0.1	●
8-204-32-9	10-100	0.1	●
8-206-32-9	20-200	0.2	●
8-203-32-9	30-300	0.2	●

12 Channel

Cat. No.	Volume Range (µL)	Increments (µL)	Colour Code
8-401-32-9	0.5-10	0.02	●
8-405-32-9	2-20	0.02	●
8-402-32-9	5-50	0.1	●
8-404-32-9	10-100	0.1	●
8-406-32-9	20-200	0.2	●
8-403-32-9	30-300	0.2	●



We also offer
the compatible
AHN myTip Effortless
EL³ tips and filter tips
for our AHN myPette EL³
pipettes



© 2024 AHN Biotechnologie GmbH. Version: 2408