

ahn myTip®

Tecan® Type Robotic Tips



Exceptional results when paired with LiHa® Arms.

Versatility - whether you are working with the Freedom EVO®100, 150, or 200 with Liquid Handling Arm, or the Freedom EVO®75 8+1, when it comes to automated liquid handling, nothing surpasses the AHN Robotic Tips for Tecan®.

Premium-grade virgin polypropylene - durability and reliability

Tips options - carbon-embedded conductive black or clear tip

Polyethylene filter - effectively prevents aerosol contamination

Made from USP class VI resins - and certified free of Endotoxin, Heavy Metals, RNase, DNase and DNA free

Pick from non-sterile, sterile and filtered sterile formats - depending on your application

Manufactured in ISO-certified 9001 and 13485, 100K cleanroom facility - to meet the highest quality standards

Validated - on the actual instrument for which the tip is created, to ensure seamless integration

Colour-coded polycarbonate tray - quickly identify volumes

Packaging - blister packs / eco-friendly flip card box















Ordering information

Tecan® type Robotic Tips






















Volume	Length	Colour Tip	Colour Code	Conductive	Sterile	Pack Type	Sales Unit	Cat. No.	
myTip® Tecan® type tip, 10 µL	31.70 mm	clear				blister	case / 24x96	1-302-B24-0	
		clear				flip box	case / 5x10x2x96	1-302-F50-0	
		black	■	■		blister	case / 24x96	1-304-B24-7	
		black	■	■		flip box	case / 5x10x2x96	1-304-F50-7	
		clear			■		blister	case / 24x96	2-302-B24-0
		black			■	■	blister	case / 24x96	2-304-B24-7
myTip® Tecan® type tip, 50 µL	58.05 mm	clear				blister	case / 24x96	1-312-B24-0	
		clear				flip box	case / 5x10x2x96	1-312-F50-0	
		black	■	■		blister	case / 24x96	1-314-B24-7	
		black	■	■		flip box	case / 5x10x2x96	1-314-F50-7	
		clear				■	blister	case / 24x96	2-312-B24-0
		black			■	■	blister	case / 24x96	2-314-B24-7

TECAN® & LiHa® are registered trademarks of Tecan Trading AG



Volume	Length	Colour Tip	Colour Code	Conductive	Sterile	Pack Type	Sales Unit	Cat. No.	
myTip® Tecan® type tip, 200 µL	58.40 mm	clear				blister	case / 24x96	1-327-B24-0	
		clear				flip box	case / 5x10x2x96	1-327-F50-0	
		black				blister	case / 24x96	1-329-B24-7	
		black				flip box	case / 5x10x2x96	1-329-F50-7	
		clear					blister	case / 24x96	2-327-B24-0
		black					blister	case / 24x96	2-329-B24-7
myTip® Tecan® type tip, 1000 µL	95.80 mm	clear				blister	case / 24x96	1-402-B24-0	
		clear				flip box	case / 5x10x2x96	1-402-F50-0	
		black				blister	case / 24x96	1-404-B24-7	
		black				flip box	case / 5x10x2x96	1-404-F50-7	
		clear					blister	case / 24x96	2-402-B24-0
		black					blister	case / 24x96	2-404-B24-7

Tecan® Type Robotic Filter Tips

Volume	Length	Colour Tip	Colour Code	Conductive	Sterile	Pack Type	Sales Unit	Cat. No.	
myTip® Tecan® type tip, 10 µL	31.70 mm	clear				blister	case / 24x96	1-303-B24-0	
		clear				flip box	case / 5x10x2x96	1-303-F50-0	
		black				blister	case / 24x96	1-305-B24-7	
		black				flip box	case / 5x10x2x96	1-305-F50-7	
		clear					blister	case / 24x96	2-303-B24-0
		black					blister	case / 24x96	2-305-B24-7
myTip® Tecan® type tip, 50 µL	58.05 mm	clear				blister	case / 24x96	1-313-B24-0	
		clear				flip box	case / 5x10x2x96	1-313-F50-0	
		black				blister	case / 24x96	1-315-B24-7	
		black				flip box	case / 5x10x2x96	1-315-F50-7	
		clear					blister	case / 24x96	2-313-B24-0
		black					blister	case / 24x96	2-315-B24-7
myTip® Tecan® type tip, 200 µL	58.40 mm	clear				blister	case / 24x96	1-328-B24-0	
		clear				flip box	case / 5x10x2x96	1-328-F50-0	
		black				blister	case / 24x96	1-330-B24-7	
		black				flip box	case / 5x10x2x96	1-330-F50-7	
		clear					blister	case / 24x96	2-328-B24-0
		black					blister	case / 24x96	2-330-B24-7
myTip® Tecan® type tip, 1000 µL	95.80 mm	clear				blister	case / 24x96	1-403-B24-0	
		clear				flip box	case / 5x10x2x96	1-403-F50-0	
		black				blister	case / 24x96	1-405-B24-7	
		black				flip box	case / 5x10x2x96	1-405-F50-7	
		clear					blister	case / 24x96	2-403-B24-0
		black					blister	case / 24x96	2-405-B24-7

TECAN® & LiHa® are registered trademarks of Tecan Trading AG