

# ahn myLab<sup>®</sup>

## Single Plate Magnetic Stirrer SM-01

**Usage** - in laboratories, pharmacies and schools, mixing liquid industry

**Maximum speed** - from 15 to 1500 rpm with an interval of 50 rpm / software-driven gradual acceleration

**Time setting** - ranging from 1 - 99 minutes / continuous operation feature

**Improved mixing** - pulse mode where the mixing direction changes every 30 seconds

**Maintenance-free design** - water spill, corrosion and chemically resistant / includes no moving parts / high-quality stirring bar

**User-friendly interface and design** - digital display showing the real-time speed and the remaining time / last run memory function / extremely thin (only 16.0mm thick) / flat surface enables easy wiping and cleaning



1500 rpm



## Ordering information

Description	Cat. No.	Specifications	
AHN myLab <sup>®</sup> SM-01 Single Plate Magnetic Stirrer 1500 rpm	7-300-00-0	Max. stirring quantity	800 mL
		Recommended stir bar length	25mm
		Protection class	IP65
		Dimensions (L x W x H)	156 x 110 x 13.2mm
		Weight	400 g
Optional accessories		Cat. No.	
Magnetic stirring bar - PTFE coated, 25 x 8mm	CRS-MAG		
Delivery Package			
1 pc. Magnetic Stirrer 1500 rpm			
1 pc. Stir bar (PTFE coated, 25 x 8mm)			
1 pc. Instruction manual			
1 pc. Warranty card			

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com



www.ahn-bio.de

