

AHN myPette® SP

Serological Pipettes
1.0 mL - 25 mL



Quick volume check - glass-clear polystyrene (PS) with graduation marks

Maintain sample integrity - sterile, individually packed

Convenient identification - colour-coded cotton-wool plug, standard colour code

Accurate, minimal residual wetting - weldless one piece pipette

Additional capacity - availability of negative graduation

Display of dispensed and residual volumes - 5, 10 and 25 mL pipettes provide reverse double scales



Description	Length	Volume
	271.0 mm	1.0 mL
	271.0 mm	2.0 mL
	294.0 mm	5.0 mL
	294.0 mm	10.0 mL
	297.0 mm	25.0 mL

Ordering information

Pos. Volume	Neg. Volume	Colour Code	Graduation	Sterilized	Pack Type	Sales Unit	Cat. No.
1.0 mL	0.3 mL		0.01 mL	■	Ind. wrapped	1x2000	8-840-00-0
				■	Bag/25 pcs.	160x25	8-840-00-1
2.0 mL	0.4 mL		0.02 mL	■	Ind. wrapped	1x2000	8-840-01-0
				■	Bag/25 pcs.	160x25	8-840-01-1
5.0 mL	3.0 mL		0.10 mL	■	Ind. wrapped	1x1000	8-840-02-0
				■	Bag/25 pcs.	60x25	8-840-02-1
10.0 mL	3.0 mL		0.10 mL	■	Ind. wrapped	1x800	8-840-03-0
				■	Bag/25 pcs.	40x25	8-840-03-1
25.0 mL	2.6 mL		0.20 mL	■	Ind. wrapped	1x400	8-840-04-0
				■	Bag/10 pcs.	40x10	8-840-04-1

