

COMPATIBILITY OF

AHN myTip FT
Standard Filter Tips

AHN myTip FT Standard Filter Tips Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN											
Brand		pipet4u [®] pro				pipet4u [®] Revolution							
Model		variable/multi channel				variable/multi channel							
		8 Channel		12 Channel		8 Channel		12 Channel					
Volume Range / Cat. No.		0.5-10 µl / 8-201-20-9	5-50 µl / 8-202-20-9	3.0-300 µl / 8-203-20-9	0.5-10 µl / 8-401-20-9	5-50 µl / 8-402-20-9	30-300 µl / 8-403-20-9	0.2-10 µl / 8-520-20-0	2.5-50 µl / 8-521-20-0	15-300 µl / 8-522-20-0	0.2-10 µl / 8-530-20-0	2.5-50 µl / 8-531-20-0	15-300 µl / 8-532-20-0
AHN myTip FT Filter Tips													
Volume	Cat.No.												
0.1-10 µl	2-001-70-0	■			■				■			■	
	2-001-C5-0	■			■				■			■	
0.1-10 µl XL	2-061-70-0	■			■				■			■	
	2-061-C5-0	■			■				■			■	
0.5-10/20 µl	2-051-70-0	■			■				■			■	
	2-051-C5-0	■			■				■			■	
1-20/30 µl	2-116-70-0	■	■	■	■	■	■		■	■	■	■	■
	2-116-C5-0	■	■	■	■	■	■		■	■	■	■	■
1-50 µl	2-111-70-0	■	■	■	■	■	■		■	■	■	■	■
	2-111-C5-0	■	■	■	■	■	■		■	■	■	■	■
1-100 µl	2-112-70-0	■	■	■	■	■	■		■	■	■	■	■
	2-112-C5-0	■	■	■	■	■	■		■	■	■	■	■
1-200 µl	2-113-70-0	■	■	■	■	■	■		■	■	■	■	■
	2-113-C5-0	■	■	■	■	■	■		■	■	■	■	■
5-300 µl	2-115-70-0	■	■	■	■	■	■		■	■	■	■	■
	2-115-C5-0	■	■	■	■	■	■		■	■	■	■	■

Filter Tips Compatibility
AHN Multi Channel Pipettes





Uthleber Weg 14
D-99734 Nordhausen
Germany

Tel.: +49(0)3631/46594-04
Fax: +49(0)3631/46594-10

E-Mail: info@ahn-bio.de
Web: www.ahn-bio.de