

COMPATIBILITY OF

AHN myTip LFT
Low Retention Filter Tips

AHN myTip LFT Low Retention Filter Tips Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN			
Brand		pipet4u® pro		pipet4u® Revolution	
Model		variable/multi channel		variable/multi channel	
		8 Channel	12 Channel	8 Channel	12 Channel
Volume Range / Cat. No.		0.5-10 µl / 8-201-20-9 5-50 µl / 8-202-20-9 30-300 µl / 8-203-20-9	0.5-10 µl / 8-401-20-9 5-50 µl / 8-402-20-9 30-300 µl / 8-403-20-9	0.2-10 µl / 8-520-20-0 2.5-50 µl / 8-521-20-0 15-300 µl / 8-522-20-0	0.2-10 µl / 8-530-20-0 2.5-50 µl / 8-531-20-0 15-300 µl / 8-532-20-0
AHN myTip LFT Filter Tips					
Volume	Cat.No.				
0.1-10 µl	5-001-C5-0	■	■	■	■
0.1-10 µl XL	5-061-C5-0	■	■	■	■
0.5-10/20 µl	5-051-C5-0	■	■	■	■
1-20/30 µl	5-116-C5-0	■	■	■	■
1-50 µl	5-111-C5-0	■	■	■	■
1-100 µl	5-112-C5-0	■	■	■	■
1-200 µl	5-113-C5-0	■	■	■	■

Filter Tips Compatibility
AHN Multi Channel Pipettes



metal-free

autoclavable
121° C | 15 min





Uthleber Weg 14
D-99734 Nordhausen
Germany

Tel.: +49(0)3631/46594-04
Fax: +49(0)3631/46594-10

E-Mail: info@ahn-bio.de
Web: www.ahn-bio.de