

COMPATIBILITY OF

AHN myTip LT
Low Retention Tips

AHN myTip Low Retention Tips Compatibility with AHN Single Channel Pipettes

Tips Compatibility
AHN Single Channel Pipettes

| Pipettor | | AHN | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|-------------------------|--|--------------------------|------------------------|----------------------|----------------------|------------------------|------------------------|------------------------|--------------------------|--------------------------|--------------------------|-------------------------------------|-------------------|--------------------|--------------------|--------------------|--------------------|---------------------|---------------------|---------------------|---------------------|----------------------|----------------------|----------------------|----------------------|-----------------------|------------------------|-----------------------|------------------------|-------------------------|--------------------------|
| Brand | | pipet4u [®] pro | | | | | | | | | | pipet4u [®] Revolution | | | | | | | | | | | | | | | | | | | |
| Model | | variable/single channel | | | | | fix/single channel | | | | | variable/ electronic single channel | | | | | | | | | | | | | | | | | | | |
| Volume Range / Cat. No. | | 0.1-2.5 µl / 8-100-20-9 | 0.5-10 µl / 8-101-20-9 | 2-20 µl / 8-102-20-9 | 5-50 µl / 8-103-20-9 | 10-100 µl / 8-104-20-9 | 20-200 µl / 8-105-20-9 | 50-200 µl / 8-105-21-9 | 100-1000 µl / 8-106-20-9 | 200-1000 µl / 8-106-21-9 | 500-5000 µl / 8-107-20-9 | 2000-10000 µl / 8-108-20-9 | 5 µl / 8-001-20-9 | 10 µl / 8-002-20-9 | 20 µl / 8-003-20-9 | 25 µl / 8-004-20-9 | 50 µl / 8-007-20-9 | 100 µl / 8-008-20-9 | 200 µl / 8-009-20-9 | 250 µl / 8-010-20-9 | 500 µl / 8-011-20-9 | 1000 µl / 8-012-20-9 | 2000 µl / 8-013-20-9 | 2500 µl / 8-014-20-9 | 5000 µl / 8-015-20-9 | 10000 µl / 8-016-20-9 | 0.2-10 µl / 8-500-20-0 | 5-100 µl / 8-501-20-0 | 10-200 µl / 8-502-20-0 | 50-1000 µl / 8-503-20-0 | 100-5000 µl / 8-505-20-0 |
| Volume | | AHN myTip LT Tips | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Cat-No. | | 4-000-81-0 | 4-000-C5-0 | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.1-10 µl | | ■ | ■ | ■ | | | | | | | | | | | | | | | | | | | | | | ■ | | | | | |
| 0.1-10 µl XL | | ■ | ■ | ■ | | | | | | | | | | | | | | | | | | | | | | ■ | | | | | |
| 0.5-10/20 µl | | | ■ | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1-200 µl | | | | | ■ | ■ | ■ | ■ | | | | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | | | |
| 100-1000 µl long | | | | | | | | ■ | ■ | | | | | | | | | ■ | ■ | | | | | | | | | | | ■ | ■ |



AHN myTip Low Retention Tips Compatibility with AHN Multi Channel Pipettes

| Pipettor | | AHN | | | | | | | | | | | |
|-------------------------|------------|------------------------|----------------------|------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|------------------------|
| Brand | | pipet4u® pro | | | | pipet4u® Revolution | | | | | | | |
| Model | | variable/multi channel | | | | variable/multi channel | | | | | | | |
| | | 8 Channel | | 12 Channel | | 8 Channel | | 12 Channel | | | | | |
| Volume Range / Cat. No. | | 0.5-10 µl / 8-201-20-9 | 5-50 µl / 8-202-20-9 | 30-300 µl / 8-203-20-9 | 0.5-10 µl / 8-401-20-9 | 5-50 µl / 8-402-20-9 | 30-300 µl / 8-403-20-9 | 0.2-10 µl / 8-520-20-0 | 2.5-50 µl / 8-521-20-0 | 15-300 µl / 8-522-20-0 | 0.2-10 µl / 8-530-20-0 | 2.5-50 µl / 8-531-20-0 | 15-300 µl / 8-532-20-0 |
| AHN myTip LT Tips | | | | | | | | | | | | | |
| Volume | Cat.No. | | | | | | | | | | | | |
| 0.1-10 µl | 4-000-81-0 | ■ | | | ■ | | | ■ | | | ■ | | |
| | 4-000-C5-0 | ■ | | | ■ | | | ■ | | | ■ | | |
| 0.1-10 µl XL | 4-060-81-0 | ■ | | | ■ | | | ■ | | | ■ | | |
| | 4-060-C5-0 | ■ | | | ■ | | | ■ | | | ■ | | |
| 0.5-10/20 µl | 4-050-81-0 | ■ | | | ■ | | | ■ | | | ■ | | |
| | 4-050-C5-0 | ■ | | | ■ | | | ■ | | | ■ | | |
| 1-200 µl | 4-110-81-0 | ■ | ■ | | ■ | ■ | | ■ | ■ | | ■ | ■ | |
| | 4-110-C5-0 | ■ | ■ | | ■ | ■ | | ■ | ■ | | ■ | ■ | |

Tips Compatibility
AHN Multi Channel Pipettes



metal-free

autoclavable
121° C | 15 min





Uthleber Weg 14
D-99734 Nordhausen
Germany

Tel.: +49(0)3631/46594-04
Fax: +49(0)3631/46594-10

E-Mail: info@ahn-bio.de
Web: www.ahn-bio.de