

COMPATIBILITY OF

AHN myTip Standard Tips

AHN myTip Macro Tips

AHN myTip Compatibility with AHN Single Channel Pipettes

Pipettor		AHN													
Brand	pipet4u® pro										pipet4u® Revolution				
Model	variable/single channel					fix/single channel					variable/electronic single channel				
Volume Range / Cat. No.	0.1-2.5 µl / 8-100-20-9 0.5-10 µl / 8-101-20-9 2-20 µl / 8-102-20-9 5-50 µl / 8-103-20-9 10-100 µl / 8-104-20-9 20-200 µl / 8-105-20-9 50-200 µl / 8-105-21-9 100-1000 µl / 8-106-20-9 200-1000 µl / 8-106-21-9 500-5000 µl / 8-107-20-9 2000-10000 µl / 8-108-20-9					5 µl / 8-001-20-9 10 µl / 8-002-20-9 20 µl / 8-003-20-9 25 µl / 8-004-20-9 50 µl / 8-007-20-9 100 µl / 8-008-20-9 200 µl / 8-009-20-9 250 µl / 8-010-20-9 500 µl / 8-011-20-9 1000 µl / 8-012-20-9 2000 µl / 8-013-20-9 2500 µl / 8-014-20-9 5000 µl / 8-015-20-9 10000 µl / 8-016-20-9					0.2-10 µl / 8-500-20-0 5-100 µl / 8-501-20-0 10-200 µl / 8-502-20-0 50-1000 µl / 8-503-20-0 100-5000 µl / 8-505-20-0				
AHN myTip Standard Tips															
Volume	Cat.No.														
0.1-10 µl	1-000-81-0 1-000-C5-0 1-001-C5-0	■	■	■	■								■	■	
0.1-10 µl XL	1-060-81-0 1-060-C5-0 1-061-C5-0	■	■	■	■								■	■	
0.5-10/20 µl	1-050-81-0 1-050-C5-0 1-051-C5-0		■	■											
1-200 µl	1-110-60-0 1-110-60-1 1-110-81-0 1-110-81-1 1-110-C5-0 1-110-C5-1 1-111-C5-0			■	■	■	■	■	■	■	■		■	■	
1-200 µl Gelfoader	1-130-70-0 1-130-C5-0 1-131-C5-0			■	■	■	■	■	■	■			■	■	
5-300 µl	1-103-81-0 1-103-C5-0 1-104-C5-0			■	■	■	■	■	■	■			■	■	
100-1000 µl short	1-209-80-2 1-209-C5-2 1-207-C5-2							■	■				■	■	
100-1000 µl long	1-200-80-0 1-200-80-2 1-200-C5-0 1-200-C5-2 1-201-C5-0									■	■	■	■	■	
AHN myTip Macro Tips															
1-5 ml	1-210-84-0 1-214-84-0														
1-5 ml	1-211-84-0 1-215-84-0													■	
1-5.5 ml	1-212-84-0 1-216-84-0														
1-10 ml	1-213-50-0 1-217-50-0												■	■	



metal-free

autoclavable
121° C | 15 min



AHN myTip Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN			
Brand		pipet4u® pro		pipet4u® Revolution	
Model	Volume Range / Cat. No.	variable/multi channel		variable/multi channel	
		8 Channel	12 Channel	8 Channel	12 Channel
		0.5-10 µl / 8-201-20-9	0.5-10 µl / 8-401-20-9	0.2-10 µl / 8-520-20-0	0.2-10 µl / 8-530-20-0
		5-50 µl / 8-202-20-9	5-50 µl / 8-402-20-9	2.5-50 µl / 8-521-20-0	2.5-50 µl / 8-531-20-0
		30-300 µl / 8-203-20-9	30-300 µl / 8-403-20-9	15-300 µl / 8-522-20-0	15-300 µl / 8-532-20-0
AHN myTip Standard Tips					
Volume	Cat.No.				
0.1-10 µl	1-000-81-0	■	■	■	■
	1-000-C5-0	■	■	■	■
	1-001-C5-0	■	■	■	■
0.1-10 µl XL	1-060-81-0	■	■	■	■
	1-060-C5-0	■	■	■	■
	1-061-C5-0	■	■	■	■
0.5-10/20 µl	1-050-81-0	■	■	■	■
	1-050-C5-0	■	■	■	■
	1-051-C5-0	■	■	■	■
1-200 µl	1-110-60-0	■	■	■	■
	1-110-60-1	■	■	■	■
	1-110-81-0	■	■	■	■
	1-110-81-1	■	■	■	■
	1-110-C5-0	■	■	■	■
	1-110-C5-1	■	■	■	■
	1-111-C5-0	■	■	■	■
1-200 µl Gelfoader	1-130-70-0	■	■	■	■
	1-130-C5-0	■	■	■	■
	1-131-C5-0	■	■	■	■
5-300 µl	1-103-81-0	■	■	■	■
	1-103-C5-0	■	■	■	■
	1-104-C5-0	■	■	■	■

Tips Compatibility
AHN Multi Channel Pipettes



metal-free

autoclavable
121° C | 15 min





Uthleber Weg 14
D-99734 Nordhausen
Germany

Tel.: +49(0)3631/46594-04
Fax: +49(0)3631/46594-10

E-Mail: info@ahn-bio.de
Web: www.ahn-bio.de