

ahn myLab®

Clinical Centrifuge CLC-01

Usage - designed for blood centrifugation and other clinical applications

Multiple rotor and tube options - 6x10 mL (swing out rotor) / 8x15 mL (fixed angle rotor) / 16x10 mL (fixed angle rotor) / with each rotor, the appropriate tube adaptors, sleeves, spacers are included - 10 mL, 9 mL, 8 mL, 6 mL, 4 mL, 3 mL and 2 mL

Maximum Speed - up to 4000 rpm / 2270 x g varies per rotor

Precise operation - speed intervals of 10 rpm, starting from 500 rpm

Memory feature - can save up to 99 user-defined programs (protocols)

Increased user safety and comfort - automatic safety brake system preventing the lid from opening upon operation / imbalance detection

Robust design - brushless DC motor enables maintenance-free and long-lasting operation, perfect for extended runs

Comfortable use - last run memory function / user-friendly interface for speed and time setting / quick spin feature



4000 rpm

Ordering information

Specifications

Rotor Type: Swing out		Rotor Type: Fixed angle		Rotor Type: Fixed angle	
Capacity: 6x10 mL		Capacity: 8x15 mL		Capacity: 16x10 mL	
Max RCF: 2270 x g		Max RCF: 1950 x g		Max RCF: 2180 x g	
Cat. No.: CRC-4-6S		Cat. No.: CRC-4-8		Cat. No.: CRC-4-16	

Description	Cat. No.
AHN myLab® CLC-01 Clinical Centrifuge 4000 rpm	7-010-00-0

Delivery Package
1 pc. Clinical centrifuge 4000 rpm
1 pc. Power supply adaptor
1 pc. T-Allen Key
1 pc. Swing out rotor
2 x 6 pcs. Tube holders
1 pc. Instruction manual
1 pc. Warranty card

Specifications	
Imbalance detection	yes, auto cut off
Timer setting	1-999 mins / infinite mode
Weight	7.7 kg (without rotor)
Dimensions (W x D x H)	355x415x173mm

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com

www.ahn-bio.de

