

AHN myLab® TR-01

Tube Roller 30 rpm



Usage - mixing of blood samples and viscous substances

Capacity - 9 x 5 mL tubes


Fixed speed - 30 rpm with an accuracy of +/-2 rpm

Optimal mixing - gentle rocking and rolling movements at a tilt angle of ± 3 degree

Easy maintenance and adjustment - tool free removal and installation of rollers for effortless cleaning and replacing / easy adjustment of spacing to use with different tube diameter

Robust and space-saving design - robust construction for stable operation / small footprint requiring minimum of bench space

Tool free removal and
installation of rollers



Ordering information

Description

AHN myLab® TR-01 Tube Roller 30 rpm

Cat. No.

7-200-00-0

Delivery package

- 1 pc. Tube Roller 30 rpm
- 1 pc. Power supply adaptor
- 1 pc. Instruction manual
- 1 pc. Warranty card

Specifications

| | |
|-----------------------|-----------------------|
| Run time | Continuous / infinite |
| Noise level | Noiseless operation |
| Roller size | 30 x 330 mm |
| Max. loading capacity | 1.5 kg |
| Protection class | IP 21 |
| Dimensions (WxDxH) | 475 x 201 x 80 mm |
| Weight | 2.0 kg |

