

# ahn myTube CT



Centrifuge Tubes 15 mL | 50 mL

**Elimination of contamination of sample** - high quality raw material meets purity standards

**Centrifuge Tube with screw flat cap**

**Easy sample identification** - large writing area on tube

**Quick and simple volume check** - precise graduation marks

**No loss of sample** - tube resists high centrifugal force and temperature variations (-80°C to +121°C)

**Centrifuge specifications:**

up to 17.000 RCF - 15 mL centrifuge tube

up to 20.000 RCF - 50 mL centrifuge tube

up to 7.000 RCF - 50 mL centrifuge tube (self-standing)

**Comfortable handling** - the 3/4th turn closure enables one-handed operation by opening and tightening

**Maintain sample integrity** - sterile version available (by gamma irradiation)

**Sterile handling** - centrifuge tubes are autoclavable (opened 121°C / 15 min)

**Versions** - 50 mL tubes are available as non-skirted or self-standing (skirted)



metal free

autoclavable  
121°C  
15 min

## Ordering information

Volume	Height	Colour	Sterilised	Pack Type	Sales Unit	Cat. No.
15 mL	120.0 mm	clear/grey		bag	20x25	3-254-25-8
		clear/grey		rack	20x25	3-255-25-8
		clear/grey	■	bag	20x25	3-256-25-8
		clear/grey	■	rack	20x25	3-257-25-8
50 mL	115.0 mm	clear/grey		bag	20x25	3-270-25-8
		clear/grey		rack	20x25	3-271-25-8
		clear/grey	■	bag	20x25	3-272-25-8
		clear/grey	■	rack	20x25	3-273-25-8
50 mL self-standing	115.0 mm	clear/grey		bag	20x25	3-280-25-8
		clear/grey	■	bag	20x25	3-281-25-8



Non-skirted



Self-standing (skirted)

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com



www.ahn-bio.de

