

AHN myLab® CLC-01

Clinical Centrifuge 4000 rpm

Accessories



Rotor sets (consisting of rotor and appropriate adaptors, sleeves, spacers)

Rotor Set 1		Rotor Set 2		Rotor Set 3	
<p>Rotor type Swing out Capacity 6x10 mL Max RCF 2270 x g Cat. No. 7-010-03-0</p>		<p>Rotor type Fixed angle Capacity 8x15 mL Max RCF 1950 x g Cat. No. 7-010-01-0</p>		<p>Rotor type Fixed angle Capacity 16x10 mL Max RCF 2180 x g Cat. No. 7-010-02-0</p>	
<p>Tube Adaptor 8 / 9 / 10 mL Quantity 6 pcs.</p>		<p>Tube Sleeve 6 mL Quantity 8 pcs.</p>		<p>Tube Sleeve 9 / 10 mL Quantity 16 pcs.</p>	
<p>Tube Sleeve 6 mL Quantity 6 pcs.</p>		<p>Tube Spacer 2 / 3 / 4 mL Quantity 8 pcs.</p>		<p>Tube Sleeve 6 mL Quantity 16 pcs.</p>	
<p>Tube Adaptor 2 / 3 / 4 mL Quantity 6 pcs.</p>				<p>Tube Spacer 2 / 3 / 4 mL Quantity 16 pcs.</p>	
<p>Tube Sleeve 2 / 3 / 4 mL Quantity 6 pcs.</p>					



+49(0)3631/46594-04
+49(0)3631/46594-10
info@ahn-bio.de

Uthleber Weg 14
D-99734 Nordhausen
www.ahn-bio.de

ahn®
biotechnologie GmbH