

ahn myLab VT-01

Mini Vortex Mixer 4500 rpm

Usage - quick mixing, re-suspension of samples in individual test tubes, up to a diameter of 30 mm

Maximum Speed - adjustable up to 4500 rpm

Robust design - brushless DC motor enables maintenance free and long lasting operation

Stable operation - through aluminium die cast base, suction cups and inbuilt counter balance system

Very small footprint saves valuable bench space - suitable for use in hoods and cold rooms



Soft silicone cup



Speed knob



Ordering information

Description	Cat. No.
AHN myLab® VT-01 Mini Vortex Mixer 4500 rpm	7-100-00-0

Specifications	
Max. tube diameter	30 mm
Orbital diameter	4.5 mm
Max. loading capacity	500 g
Press activation	yes
Protection class	IP 43
Dimensions (WxDxH)	100x100x67 mm
Weight	0.55 kg

Delivery Package
1 pc. Mini Vortex Mixer 4500 rpm
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card

