

# ahn myPlate PCR





PCR Plates 96 x 0.1 mL | 96 x 0.2 mL

**Non-skirted** - for optimal fit in various thermocyclers /can be easily cut into 16, 24, 32 or 48-well pieces

**Optical measurement** - high transparency of PCR plates allows the use for applications where the optical qPCR signal comes from the bottom of the plate

**Improved sample recovery** - the V-shape bottom guarantees easier sample recovery



**Low-binding** - reducing the adhesion of proteins to the very minimum

**Sterile handling** - each lot is manufactured under strict sterile conditions

**Sealable** - using optical disposable adhesive for secure closure of PCR plate

**Versions** - available as low profile PCR plates with 0.1 mL wells and as regular profile plates with 0.2 mL well capacity

**Compatibility chart** available on request

metal free

autoclavable  
121°C  
15 min

## Ordering information

Volume	Height	Colour	Pack Type	Sales Unit	Cat. No.
96 x 0.1 mL	16.6 mm	clear	bag	8x25	3-046-00-0
96 x 0.2 mL	21.8 mm	clear	bag	8x25	3-045-00-0

Accessory	Pack Type	Sales Unit	Cat. No.
Optical disposable adhesive	bag	1x100	3-047-00-0

