

ahn myTip FT

Pipette Filter Tips



No sample contamination - filters are made of virgin high quality hydrophobic polyethylene without additives

No pipette contamination - filters block aerosols and liquids (vaporious, radioactive, biohazardous or corrosive samples) through micropore filter technology (MPS technology, filter pore size 5-40 μm)

Universal fit and tip ejection - optimised cone geometry for a secure fit

Quick and easy volume check - graduation marks on tip at standard volumes

Precise pipetting without loss of sample - conical tip opening prevents drop formation at the end of the tip


Elimination of sample contamination - tips made of autoclavable, non-wettable, virgin polypropylene

Easy identification of pipette tip - colour coded racks allow a quick and easy determination of the pipette tip and the suitable pipette

Maintain sample integrity - sterile version available (beta irradiation)

Versions - available in bags and rack boxes









Compatibility list available on request



Filter	Average pore size	Material	Characteristic	Used in following filter tip
	18-40 μm	PE	hydrophobic	10 μL 10 μL XL
	5-20 μm	PE	hydrophobic	20 μL
	5-20 μm	PE	hydrophobic	100 μL
	5-20 μm	PE	hydrophobic	200 μL 300 μL
	18-40 μm	PE	hydrophobic	1000 μL 1250 μL



Ordering information

Volume	Length	Colour Tip	Colour Code	Sterilised	Pack Type	Sales Unit	Cat. No.
 AHN myTip® FT 10 µL	31.20 mm	clear			bag	case / 20 x 1000 tips	2-002-81-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	2-003-C5-0
 AHN myTip® FT 10 µL XL	45.80 mm	clear			bag	case / 20 x 1000 tips	2-062-81-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	2-063-C5-0
 AHN myTip® FT 20 µL	53.40 mm	clear			bag	case / 10 x 1000 tips	2-117-80-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	2-118-C5-0
 AHN myTip® FT 100 µL	53.40 mm	clear			bag	case / 10 x 1000 tips	2-119-80-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	2-120-C5-0
 AHN myTip® FT 200 µL	53.40 mm	clear			bag	case / 10 x 1000 tips	2-127-80-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	2-128-C5-0
 AHN myTip® FT 300 µL	59.35 mm	clear			bag	case / 10 x 1000 tips	2-134-80-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	2-135-C5-0
 AHN myTip® FT 1000 µL	84.30 mm	clear			bag	case / 10 x 1000 tips	2-202-80-0
		clear	■	■	racked	case / 4 x 8 x 96 tips	2-203-C4-0
 AHN myTip® FT 1250 µL	97.60 mm	clear			bag	case / 10 x 1000 tips	2-204-80-0
		clear	■	■	racked	case / 4 x 8 x 96 tips	2-205-C4-0

