

ahn *myTip*

Pipette Tips | Macro Tips



ahn *myTip* Pipette Tips

Optimised fit and easy tip ejection - improved cone geometry for a secure fit

Quick and easy volume check - graduation marks on tip at standard volumes

Precise pipetting without loss of sample - conical tip opening prevents drop formation at the end of the tip


High quality material - tips made of autoclavable, non-wettable, virgin polypropylene

Easy identification of pipette tip - colour coded racks allow a quick and easy determination of the pipette tip and the suitable pipette

Maintain sample integrity - sterile version available (beta irradiation)

Versions - available in bags and rack boxes

Compatibility list available on request



metal
free



ahn *myTip* MaT Macro Tips

Optimised cone geometry - ensuring perfect fit on the defined pipette

Precision without loss of sample - conical tip opening to avoid drop formation

Elimination of sample contamination - tips made of autoclavable, non-wettable, virgin polypropylene material

Maintain sample integrity - sterile version available (beta irradiation)











Packaging - available in bags




Compatibility list available on request

metal
free



Ordering information

Volume	Length	Colour Tip	Colour Code	Sterilised	Pack Type	Sales Unit	Cat. No.
 AHN myTip® 10 µL	31.20 mm	clear			bag	case / 20 x 1000 tips	1-002-81-0
		clear	■		racked	case / 5 x 10 x 96 tips	1-002-C5-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	1-003-C5-0
 AHN myTip® 10 µL, XL	45.80 mm	clear			bag	case / 20 x 1000 tips	1-062-81-0
		clear	■		racked	case / 5 x 10 x 96 tips	1-062-C5-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	1-063-C5-0
 AHN myTip® 200 µL	53.40 mm	clear			bag	case / 30 x 1000 tips	1-121-86-0
		clear	■		racked	case / 5 x 10 x 96 tips	1-121-C5-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	1-122-C5-0
 AHN myTip® 200 µL	53.40 mm	yellow			bag	case / 30 x 1000 tips	1-121-86-1
		yellow	■		racked	case / 5 x 10 x 96 tips	1-121-C5-1
		yellow	■	■	racked	case / 5 x 10 x 96 tips	1-122-C5-1
 AHN myTip® Gelloader 200 µL	68.50 mm	clear			bag	case / 5 x 1000 tips	1-130-70-0
		clear	■		racked	case / 4 x 8 x 96 tips	1-130-C4-0
		clear	■	■	racked	case / 4 x 8 x 96 tips	1-131-C4-0
 AHN myTip® 300 µL	59.35 mm	clear			bag	case / 10 x 1000 tips	1-127-80-0
		clear	■		racked	case / 5 x 10 x 96 tips	1-127-C5-0
		clear	■	■	racked	case / 5 x 10 x 96 tips	1-128-C5-0
 AHN myTip® 1000 µL, short	70.10 mm	clear blue			bag	case / 10 x 1000 tips	1-209-80-2
 AHN myTip® 1000 µL	84.30 mm	clear			bag	case / 10 x 1000 tips	1-202-80-0
		clear	■		racked	case / 4 x 8 x 96 tips	1-202-C4-0
		clear	■	■	racked	case / 4 x 8 x 96 tips	1-203-C4-0
 AHN myTip® 1000 µL	84.30 mm	blue			bag	case / 10 x 1000 tips	1-202-80-2
		blue	■		racked	case / 4 x 8 x 96 tips	1-202-C4-2
		blue	■	■	racked	case / 4 x 8 x 96 tips	1-203-C4-2
 AHN myTip® 1250 µL	97.60 mm	clear			bag	case / 10 x 1000 tips	1-204-80-0
		clear	■		racked	case / 4 x 8 x 96 tips	1-204-C4-0
		clear	■	■	racked	case / 4 x 8 x 96 tips	1-205-C4-0

Volume	Length	Colour	Sterilised	Pack Type	Sales Unit	Cat. No.
 AHN myTip® MaT 5 mL	124.00 mm	clear		bag	case / 10 x 250 tips	1-218-61-0
		clear	■	bag	case / 10 x 250 tips	1-222-61-0
 AHN myTip® MaT 5.5 mL	147.50 mm	clear		bag	case / 8 x 250 tips	1-212-84-0
		clear	■	bag	case / 8 x 250 tips	1-216-84-0
 AHN myTip® MaT 10 mL	181.00 mm	clear		bag	case / 5 x 200 tips	1-220-50-0
		clear	■	bag	case / 5 x 200 tips	1-224-50-0

