

ahn myLab SM-01

Single Plate Magnetic Stirrer 1500 rpm

Usage - in laboratories, pharmacies and schools, mixing liquid industry


Maximum speed - from 15 to 1500 rpm with an interval of 50 rpm / Software driven gradual acceleration

Time setting - ranging from 1- 99 minutes / continuous operation feature

Improved mixing - pulse mode where the mixing direction changes every 30 seconds

Maintenance free design - water spill, corrosion and chemically resistant / includes no moving parts / high quality stirring bar

User-friendly interface and design - digital display showing the real-time speed and the remaining time / last run memory function / extremely thin (only 16.0 mm thick)/flat surface enables easy wiping and cleaning



Ordering information

Description		Cat. No.	Specifications	
AHN myLab® SM-01 Single Plate Magnetic Stirrer 1500 rpm		7-300-00-0	Max. stirring quantity	800 mL
			Recommended stir bar length	25 mm
Optional accessories		Cat. No.	Protection class	IP65
Magnetic stirring bar - PTFE coated, 25 x 8 mm		CRS-MAG	Dimensions (WxDxH)	20x164x16 mm
			Weight	400 g
Delivery Package				
1 pc. Magnetic Stirrer 1500 rpm				
1 pc. Stir bar (PTFE coated, 25 x 8 mm)				
1 pc. Instruction manual				
1 pc. Warranty card				

