

Thermocycler Compatibility AHN myTube PCR

PCR Strips 8 x 0.2 mL



| | |
|---|--|
| ■ | compatible |
| × | does not fit 0.5 mL cyclers |
| ○ | single tubes or strips are not recommended |
| ▲ | tubes and strips do not fit |

| Manufacturer - Applied Biosystems | | |
|-------------------------------------|-----------------------------|---------------|
| Model | Block | Compatibility |
| 310 | 96 x 0.2mL | ■ |
| 3100 | 96 x 0.2mL | ■ |
| 3130 | 96 x 0.2mL | ■ |
| 3500 Genetic Analyzer | 8-capillary | ■ |
| 3500xL Genetic Analyzer | 24-capillary | ■ |
| 3700 | 96 x 0.2mL | ■ |
| 3730 | 96 x 0.2mL | ■ |
| 7500 Real-Time PCR System | 96 x 0.2mL Standard | ■ |
| 7900HT Real-Time PCR System | 96 x 0.2mL Standard | ▲ |
| GeneAmp PCR System 2400 | 24 x 0.2mL | ■ |
| GeneAmp PCR System 2700 | 96 x 0.2mL | ■ |
| GeneAmp PCR System 2720 | 96 x 0.2mL | ■ |
| GeneAmp PCR System 5700 | 96 x 0.2mL | ■ |
| GeneAmp PCR System 9600 | 96 x 0.2mL | ■ |
| GeneAmp PCR System 9700 | 96 x 0.2mL | ■ |
| GeneAmp PCR System 9700 | Dual 96-Well 96 x 0.2mL | ■ |
| GeneAmp PCR System 9700 Aluminum | 96 x 0.2mL | ■ |
| GeneAmp PCR System 9700 Gold | 96 x 0.2mL | ■ |
| ProFlex™ PCR System | 32 x 0.2mL | ■ |
| ProFlex™ PCR System | 96 x 0.2mL | ■ |
| QuantStudio™ 12k Flex | 96 x 0.2mL | ○ |
| QuantStudio™ 6 Flex | 96 x 0.2mL | ○ |
| QuantStudio™ 7 Flex | 96 x 0.2mL | ○ |
| SimpliAmp™ | 96 x 0.2mL | ■ |
| Solid3 | 96 x 0.2mL | ■ |
| Veriti™ 96-Well Thermal Cycler | 96 x 0.2mL | ■ |
| Veriti™ Dx 96-Well Thermal Cycler | 96 x 0.2mL | ■ |
| ViiA™ 7 Real-Time PCR System | 96 x 0.2mL | ○ |
| Manufacturer - Bio-Rad Laboratories | | |
| Model | Block | Compatibility |
| C1000 | 96-Well Fast | ■ |
| C1000 | Dual 48/48 Fast | ■ |
| C1000 Touch | 96-Deep Well | ■ |
| DNA Engine | 48-48 Alpha unit 48 x 0.2mL | ■ |
| DNA Engine | 96v Alpha unit 96 x 0.2mL | ■ |
| DNA Engine Dyad | 48-48 Alpha unit 48 x 0.2mL | ■ |
| DNA Engine Dyad | 96v Alpha unit 96 x 0.2mL | ■ |
| DNA Engine Dyad Disciple | 48-48 Alpha unit 48 x 0.2mL | ■ |
| DNA Engine Dyad Disciple | 96v Alpha unit 96 x 0.2mL | ■ |
| DNA Engine Tetrad 2 | 48-48 Alpha unit 48 x 0.2mL | ■ |
| DNA Engine Tetrad 2 | 96v Alpha unit 96 x 0.2mL | ■ |

Thermocycler Compatibility

AHN myTube PCR

PCR Strips 8 x 0.2 mL



| | |
|---|--|
| ■ | compatible |
| × | does not fit 0.5 mL cyclers |
| ○ | single tubes or strips are not recommended |
| ▲ | tubes and strips do not fit |

| | | |
|---------------------|------------------|---|
| GeneCycler | 24 x 0.2mL | ■ |
| iCycler | 48 x 0.2mL | ■ |
| iCycler | 96 x 0.2mL | ■ |
| MiniOpticon | 48-12 48 x 0.2mL | ■ |
| MJ Mini | 48-12 48 x 0.2mL | ■ |
| MyCycler | 96 x 0.2mL | ■ |
| MyiQ Real-Time PCR | 96 x 0.2mL | ■ |
| MyiQ2 Real-Time PCR | 96 x 0.2mL | ■ |
| PTC-100 | 96 x 0.2mL | ■ |
| S1000 | 96-Well Fast | ■ |
| S1000 | Dual 48/48 Fast | ■ |
| T100 | 96 x 0.2mL | ■ |

Manufacturer - Eppendorf

| Model | Block | Compatibility |
|---------------------------------|------------------|---------------|
| Mastercycler | 96 x 0.2mL | ■ |
| Mastercycler | 96-77 96 x 0.2mL | ■ |
| Mastercycler ep gradient | 96 x 0.2mL | ■ |
| Mastercycler ep gradient Silver | 96 x 0.2mL | ■ |
| Mastercycler ep realplex | 96 x 0.2mL | ■ |
| Mastercycler ep realplex Silver | 96 x 0.2mL | ■ |
| Mastercycler gradient | 96 x 0.2mL | ■ |
| Mastercycler gradient | 96-77 96 x 0.2mL | ■ |
| Mastercycler nexus | 96 x 0.2mL | ■ |
| Mastercycler nexus eco | 96 x 0.2mL | ■ |
| Mastercycler nexus gradient | 96 x 0.2mL | ■ |
| Mastercycler nexus gradient eco | 96 x 0.2mL | ■ |
| Mastercycler nexus GSX1 | 96 x 0.2mL | ■ |
| Mastercycler nexus GSX1e | 96 x 0.2mL | ■ |
| Mastercycler nexus SX1 | 96 x 0.2mL | ■ |
| Mastercycler nexus SX1e | 96 x 0.2mL | ■ |
| Mastercycler Personal | 25-16 25 x 0.2mL | ■ |
| Mastercycler pro | 96 x 0.2mL | ■ |
| Mastercycler pro S | 96 x 0.2mL | ■ |
| Thermomixer comfort | 96 x 0.2mL | ■ |

Manufacturer - MJ Research Inc.

| Model | Block | Compatibility |
|-------------------------|--------------------------------|---------------|
| PTC-100 | 96 x 0.2mL | ■ |
| PTC-100 | 96 x 0.2mL AgV | ■ |
| PTC-200 DNA Engine | 30-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-200 DNA Engine | 48-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-200 DNA Engine | 96v Alpha unit 96 x 0.2mL | ■ |
| PTC-200 DNA Engine | 96v Moto Alpha Unit 96 x 0.2mL | ■ |
| PTC-220 DNA Engine Dyad | 30-48 Alpha unit 48 x 0.2mL | ■ |

Thermocycler Compatibility AHN myTube PCR

PCR Strips 8 x 0.2 mL



| | |
|---|--|
| ■ | compatible |
| × | does not fit 0.5 mL cyclers |
| ○ | single tubes or strips are not recommended |
| ▲ | tubes and strips do not fit |

| | | |
|--|--------------------------------|----------------------|
| PTC-220 DNA Engine Dyad | 48-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-220 DNA Engine Dyad | 96v Alpha unit 96 x 0.2mL | ■ |
| PTC-220 DNA Engine Dyad | 96v Moto Alpha Unit 96 x 0.2mL | ■ |
| PTC-221 DNA Engine Disciple | 30-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-221 DNA Engine Disciple | 48-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-221 DNA Engine Disciple | 96v Alpha unit 96 x 0.2mL | ■ |
| PTC-221 DNA Engine Disciple | 96v Moto Alpha Unit 96 x 0.2mL | ■ |
| PTC-225 DNA Engine Tetrad | 30-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-225 DNA Engine Tetrad | 48-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-225 DNA Engine Tetrad | 96v Alpha unit 96 x 0.2mL | ■ |
| PTC-225 DNA Engine Tetrad | 96v Moto Alpha Unit 96 x 0.2mL | ■ |
| PTC-240 DNA Engine Tetrad 2 | 30-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-240 DNA Engine Tetrad 2 | 48-48 Alpha unit 48 x 0.2mL | ■ |
| PTC-240 DNA Engine Tetrad 2 | 96v Alpha unit 96 x 0.2mL | ■ |
| PTC-240 DNA Engine Tetrad 2 | 96v Moto Alpha Unit 96 x 0.2mL | ■ |
| Manufacturer - Stratagene | | |
| Model | Block | Compatibility |
| Robocycler Gradient 40 | 96 x 0.2mL | ■ |
| Manufacturer - Whatman Biometra | | |
| Model | Block | Compatibility |
| T1 Thermocycler | 96 x 0.2mL | ■ |
| T1 Thermocycler Combi | 96-77 96 x 0.2mL | ■ |
| T3 Thermocycler | 48 x 0.2mL | ■ |
| T3 Thermocycler Combi | 48-18 48 x 0.2mL | ■ |
| T3000 Thermocycler | 48 x 0.2mL | ■ |
| T3000 Thermocycler Combi | 48-18 48 x 0.2mL | ■ |
| TGradient | 96 x 0.2mL | ■ |
| TPersonal | 48 x 0.2mL | ■ |
| TPersonal Combi | 48-18 48 x 0.2mL | ■ |
| TProfessional | 96 x 0.2mL | ■ |
| TProfessional Basic | 96 x 0.2mL | ■ |
| TProfessional Basic Gradient | 96 x 0.2mL | ■ |
| TProfessional TRIO | 48 x 0.2mL | ■ |
| TProfessional TRIO combi | 48-30 30 x 0.5mL | × |
| TProfessional TRIO combi | 48-30 48 x 0.2mL | ■ |
| TRobot | 96 x 0.2mL | ■ |
| UNO II | 96 x 0.2mL | ■ |
| UNO-thermoblock | 96 x 0.2mL | ■ |