

# ahn myLab TR-01

Tube Roller 30 rpm



**Usage** - mixing of blood samples and viscous substances

**Capacity** - 9 x 5 mL tubes

**Fixed speed** - 30 rpm with an accuracy of  $\pm 2$  rpm

**Optimal mixing** - gentle rocking and rolling movements at a tilt angle of  $\pm 3$  degree

**Easy maintenance and adjustment** - tool free removal and installation of rollers for effortless cleaning and replacing / easy adjustment of spacing to use with different tube diameter

**Robust and space-saving design** - robust construction for stable operation / small footprint requiring minimum of bench space



Tool-free removal and installation of rollers

## Ordering information

Description	Cat. No.
AHN myLab® TR-01 Tube Roller 30 rpm	7-200-00-0

Specifications	
Run time	Continuous / infinite
Noise level	Noiseless operation
Roller size	30 x 330 mm
Max. loading capacity	1.5 kg
Protection class	IP 21
Dimensions (WxDxH)	475x201x80 mm
Weight	2.0 kg

Delivery Package
1 pc. Tube Roller 30 rpm
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card

