

# ahn myDpette

## Bottle Top Dispenser

**Unique piston mechanism** - cleaning of piston and cylinder without affecting the calibration

**Dual graduation marks** - readability of the volume from both sides

**Springless valve design** - leads to smooth functioning and high chemical resistance

**Eco-friendly** - while discharging the tubes in the re-circulation mode, in fast priming, the reagent is being redirected back into the bottle, eliminating any waste of precious reagents, ensuring safety and savings

**Easy handling** - discharge tube with 360° rotation, makes the bottle label visible all the time / telescoping filling tube compatible with different sizes of reagent bottles / ergonomic grip

**Easy to calibrate** - complies with ISO 8655-5 standard




## Ordering information

Description	Sales Unit	Cat. No.
AHN myDpette Bottle Top Dispenser <b>0.25-2.5 mL</b>	Carton/1 pc.	8-700-00-0
AHN myDpette Bottle Top Dispenser <b>0.5-5.0 mL</b>	Carton/1 pc.	8-701-00-0
AHN myDpette Bottle Top Dispenser <b>1.0-10.0 mL</b>	Carton/1 pc.	8-702-00-0
AHN myDpette Bottle Top Dispenser <b>2.5-25.0 mL</b>	Carton/1 pc.	8-703-00-0
AHN myDpette Bottle Top Dispenser <b>5.0-50.0 mL</b>	Carton/1 pc.	8-704-00-0

**2**  
YEARS  
WARRANTY

autoclavable  
120°C  
20 min

Delivery Package
1 pc. Bottle top dispenser
4 pcs. Adapters (28, 38, 40 and 45 mm screw cap)
1 pc. Telescopic filling tube
1 pc. Re-circulation tube
1 pc. Calibration tool
1 pc. User manual
1 pc. Certificate of conformity
1 pc. Calibration report
1 pc. Warranty card



Double scale



Recirculation valve

