

ahn myLab CLC-01












Clinical Centrifuge 4000 rpm



Accessories

Rotor sets

(consisting of rotor and appropriate adaptors, sleeves, spacers)

Rotor Set 1		Rotor Set 2		Rotor Set 3	
Rotor Type: Swing out Capacity: 6x10 mL Max RCF: 2270 x g Cat. No.: CRC-4-6S		Rotor Type: Fixed angle Capacity: 8x15 mL Max RCF: 1950 x g Cat. No.: CRC-4-8		Rotor Type: Fixed angle Capacity: 16x10 mL Max RCF: 2180 x g Cat. No.: CRC-4-16	
Tube Adaptor: 8 / 9 / 10 mL Quantity: 6 pieces		Tube Sleeve: 6 mL Quantity: 8 pieces		Tube Sleeve: 9 / 10 mL Quantity: 16 pieces	
Tube Sleeve: 6 mL Quantity: 6 pieces		Tube Spacer: 2 / 3 / 4 mL Quantity: 8 pieces		Tube Adaptor: 6 mL Quantity: 16 pieces	
Tube Adaptor: 2 / 3 / 4 mL Quantity: 6 pieces		Tube Spacer: 2 / 3 / 4 mL Quantity: 16 pieces			
Tube Adaptor: 2 / 3 / 4 mL Quantity: 6 pieces	