

# ahn myLab VT-01

## Mini Vortex Mixer 4500 rpm

**Usage** - quick mixing, re-suspension of samples in individual test tubes, up to a diameter of 30 mm

**Maximum Speed** - adjustable up to 4500 rpm

**Robust design** - brushless DC motor enables maintenance free and long lasting operation

**Stable operation** - through aluminium die cast base, suction cups and inbuilt counter balance system

**Very small footprint saves valuable bench space** - suitable for use in hoods and cold rooms



Soft silicone cup



Speed knob



## Ordering information

Description	Cat. No.
AHN myLab® VT-01 Mini Vortex Mixer 4500 rpm	7-100-00-0

Specifications	
Max. tube diameter	30 mm
Orbital diameter	4.5 mm
Max. loading capacity	500 g
Press activation	yes
Protection class	IP 43
Dimensions (WxDxH)	100x100x67 mm
Weight	0.55 kg

Delivery Package
1 pc. Mini Vortex Mixer 4500 rpm
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card

