

ahn myPlate PCR





PCR Plates 96 x 0.1 mL | 96 x 0.2 mL

Non-skirted - for optimal fit in various thermocyclers /can be easily cut into 16, 24, 32 or 48-well pieces

Optical measurement - high transparency of PCR plates allows the use for applications where the optical qPCR signal comes from the bottom of the plate

Improved sample recovery - the V-shape bottom guarantees easier sample recovery



Low-binding - reducing the adhesion of proteins to the very minimum

Sterile handling - each lot is manufactured under strict sterile conditions

Sealable - using optical disposable adhesive for secure closure of PCR plate

Versions - available as low profile PCR plates with 0.1 mL wells and as regular profile plates with 0.2 mL well capacity


Compatibility chart available on request

metal free

autoclavable
121°C
15 min

Ordering information

Volume		Height	Colour	Pack Type	Sales Unit	Cat. No.
96 x 0.1 mL		16.6 mm	clear	bag	8x25	3-046-00-0
96 x 0.2 mL		21.8 mm	clear	bag	8x25	3-045-00-0

Accessory		Pack Type	Sales Unit	Cat. No.
Optical disposable adhesive		bag	1x100	3-047-00-0

