

Compatibility Chart

Description	Cat. No.	AHN		Gilson		Eppendorf		Brand		Thermo		Sartorius		HTL	
		Pipetku Pro / Pipetku Revolution single and multi channel pipettes	Pipetman / Pipetman G single and multi channel pipettes	Research Plus single and multi channel pipettes	Transferpette, Transferpette S single and multi channel pipettes	Finnpipette F1 / Novus manual single and multi channel pipettes and electronic single and multi channel pipettes	Proline Plus / Picus manual single and multi channel pipettes	Labmate Pro (same as Labnet BioPette) manual single and multi channel pipettes							
AHN myTip® 10 µL, clear	1-002-81-0														
AHN myTip® 10 µL, clear	1-002-C5-0														
AHN myTip® 10 µL, clear, sterile	1-003-C5-0														
AHN myTip® 10/20 µL, clear	1-062-81-0														
AHN myTip® 10/20 µL, clear	1-062-C5-0														
AHN myTip® 10/20 µL, clear, sterile	1-063-C5-0														
AHN myTip® 200 µL, clear	1-121-86-0														
AHN myTip® 200 µL, clear	1-121-C5-0														
AHN myTip® 200 µL, clear, sterile	1-122-C5-0														
AHN myTip® 200 µL, yellow	1-121-86-1														
AHN myTip® 200 µL, yellow	1-121-C5-1														
AHN myTip® 200 µL, yellow, sterile	1-122-C5-1														
AHN myTip® 200 µL, Celloclader, clear	1-130-70-0														
AHN myTip® 200 µL, Celloclader, clear	1-130-C4-0														
AHN myTip® 200 µL, Celloclader, clear, sterile	1-131-C4-0														
AHN myTip® 300 µL, clear	1-127-80-0														
AHN myTip® 300 µL, clear	1-127-C5-0														
AHN myTip® 300 µL, clear, sterile	1-128-C5-0														
AHN myTip® 1000 µL, short, clear blue	1-209-80-2														
AHN myTip® 1000 µL, clear	1-202-80-0														
AHN myTip® 1000 µL, clear	1-202-C4-0														
AHN myTip® 1000 µL, clear, sterile	1-203-C4-0														
AHN myTip® 1000 µL, blue	1-202-80-2														
AHN myTip® 1000 µL, blue, sterile	1-203-C4-2														
AHN myTip® 1250 µL, clear	1-204-80-0														
AHN myTip® 1250 µL, clear	1-204-C4-0														
AHN myTip® 1250 µL, clear, sterile	1-205-C4-0														
AHN myTip® Mat Macro Pipette Tips															
AHN myTip® Mat 5 mL, clear	1-218-61-0														
AHN myTip® Mat 5 mL, clear, sterile	1-222-61-0														
AHN myTip® Mat 5 mL, clear	1-210-84-0														
AHN myTip® Mat 5.5 mL, clear	1-212-84-0														
AHN myTip® Mat 5.5 mL, clear, sterile	1-216-84-0														
AHN myTip® Mat 10 mL, clear	1-213-50-0														
AHN myTip® Mat 10 mL, clear	1-220-50-0														
AHN myTip® Mat 10 mL, clear, sterile	1-224-50-0														

■ compatible x not compatible

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested. For additional compatibility charts, please contact info@cappahn.com.



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com



www.ahn-bio.de



Compatibility Chart

Description	Cat. No.	AHN			Gilson			Eppendorf			Brand			Thermo			Sartorius			HTL		
		Pipetku Pro / Pipetku Revolution single and multi channel pipettes	Pipetman / Pipetman G single and multi channel pipettes	Research Plus single and multi channel pipettes	Transferpette, Transferpette S single and multi channel pipettes	Finnpipette F1 / Novus manual single and multi channel pipettes and electronic single and multi channel pipettes	Proline Plus / Picus manual single channel and electronic single and multi channel pipettes	Labmate Pro (same as Labnet BioPette) manual single and multi channel pipettes														
AHN myTip® FT 10 µL, clear	2-002-81-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 10 µL, clear, sterile	2-003-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 10/20 µL, clear	2-062-81-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 10/20 µL, clear, sterile	2-063-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 20 µL, clear	2-117-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 20 µL, clear, sterile	2-118-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 100 µL, clear	2-119-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 100 µL, clear, sterile	2-120-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 200 µL, clear	2-127-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 200 µL, clear, sterile	2-128-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 300 µL, clear	2-134-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 300 µL, clear, sterile	2-135-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 1000 µL, clear	2-202-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 1000 µL, clear, sterile	2-203-C4-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 1250 µL, clear	2-204-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® FT 1250 µL, clear, sterile	2-205-C4-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT Low Retention Pipette Tips		■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 10 µL, clear	4-002-81-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 10 µL, clear, sterile	4-003-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 10/20 µL, clear	4-062-81-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 10/20 µL, clear, sterile	4-063-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 200 µL, clear	4-121-86-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 200 µL, clear, sterile	4-122-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 300 µL, clear	4-127-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 300 µL, clear, sterile	4-128-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 1000 µL, clear	4-202-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 1000 µL, clear, sterile	4-203-C4-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 1250 µL, clear	4-204-80-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LT 1250 µL, clear, sterile	4-205-C4-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT Low Retention Filter Tips		■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 10 µL, clear, sterile	5-003-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 10/20 µL, clear, sterile	5-063-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 20 µL, clear, sterile	5-118-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 100 µL, clear, sterile	5-120-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 200 µL, clear, sterile	5-128-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 300 µL, clear, sterile	5-135-C5-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 1000 µL, clear, sterile	5-203-C4-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
AHN myTip® LFT 1250 µL, clear, sterile	5-205-C4-0	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	

■ compatible x not compatible

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested. For additional compatibility charts, please contact info@cappahn.com.



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com

www.ahn-bio.de

