

ahn myTube MCT

Micro Tubes 0.5 mL | 1.5 mL | 2.0 mL



Effortless and secure handling - ergonomic cap design with hinge orientation nibs (1.5 mL / 2.0 mL) and safe click lock, allowing an easy and secure opening and closing of the tube / bottom and cap have no sharp points, gloves will not be damaged / cap of the tube is puncturable

Quick and simple volume check - tubes are equipped with precise graduation marks / 1.5 mL tube with 0.2 mL extra space and 0.5 mL tube with 0.1 mL extra space for freezing / boiling purposes

Easy sample identification - frosted marking area for easy labeling

Optical measurement and unsurpassed clarity

Elimination of contamination and loss of sample - high quality raw material meets standards of purity and is highly resistant to mechanical strain and has an excellent chemical resistance

Centrifuge specifications: up to 20.000 x g - 0.5 mL tube up to 25.000 x g - 1.5 mL / 2.0 mL tube

Maintain sample integrity - sterile version available




Sterile handling - micro tubes are autoclavable (opened 121°C / 15 min)



metal free

autoclavable
121°C
15 min

Ordering information

Volume	Height	Colour	Sterilised	Pack Type	Sales Unit	Cat. No.
0.5 mL 	31.6 mm	clear	■	bag	5x1000	3-116-C5-0
1.5 mL 	40.9 mm	clear	■	bag	16x500	3-205-80-0
		black*	■	bag	16x500	3-205-80-7
2.0 mL 	41.7 mm	clear	■	bag	16x500	3-204-80-0
		clear	■	bag	10x500	3-208-C5-0

*protection of light sensitive samples

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com



www.ahn-bio.de

