

ahn pipet4u Oasis

Pipette Controller

Comfortable and effortless pipetting - light weight pipette design /smooth push buttons for precise pipetting control / ergonomic design for easy operation / cordless working (holds charge up to 8 hours) / quiet high-powered motor

Contamination control - replaceable membrane filter to minimise risk of contamination

Speed setting - adjustment of pipetting speed on sample viscosity with 3 selectable speeds

Sterile handling - autoclavable pipette cone (120°C / 15 min)



Ordering information


| Description | Sales Unit | Cat. No. |
|--|--------------|------------|
| pipet4u® Oasis Pipette Controller 0.1-100 mL | Carton/1 pc. | 8-300-00-9 |

| Description - Spare Parts | Sales Unit | Cat. No. |
|------------------------------|------------|------------|
| 0.20 µm PTFE membrane filter | Bag/5 pcs. | 8-300-03-0 |
| 0.45 µm PTFE membrane filter | Bag/5 pcs. | 8-300-04-0 |
| Silicone adapter | Bag/1 pc. | 8-300-01-0 |
| Pipette cone | Bag/1 pc. | 8-300-02-8 |

| Delivery Package |
|---|
| 1 pc. Pipette Controller |
| 2 pcs. PTFE Membrane filters (0.2 µm / 0.45 µm) |
| 1 pc. Charger (EU) |
| 2 pcs. NiMH batteries |
| 1 pc. Wall mount bracket |
| 1 pc. Desktop stand |
| 1 pc. Instruction manual |
| 1 pc. Warranty card |



Autoclavable tip cone



Ergonomic push buttons

