



Kompatibilitätstabelle

Kat.-Nr.	Volumen	CAPP										Sartorius		Brand		Eppendorf		Finnpipette		Gilson Pipetman		Rainin		VWR, Labnet, HTL				
		CAPPBravo, ecopipette und CAPP Aero Single Channel					CAPP electronic pipettes und CAPP Aero Multi Channel					Biohit Elite Proline® Plus mLNE		Transferpette®		Research® plus		Finnpipette Digital		Gilson Pipetman		Rainin Classic, universal shaft (non LTS)		Mechanische Pipetten				
AHN myTip® Pipettenspitzen																												
1-002-81-0	AHN myTip® 10 µL, klar	■	■																									
1-002-C5-0	AHN myTip® 10 µL, klar	■	■																									
1-003-C5-0	AHN myTip® 10 µL, klar, steril	■	■																									
1-062-81-0	AHN myTip® 10 µL, XL klar	■	■																									
1-062-C5-0	AHN myTip® 10 µL, XL klar	■	■																									
1-063-C5-0	AHN myTip® 10 µL, XL steril	■	■																									
1-121-86-0	AHN myTip® 200 µL, klar		■	■	■	■																						
1-121-C5-0	AHN myTip® 200 µL, klar		■	■	■	■																						
1-222-C5-0	AHN myTip® 200 µL, klar, steril		■	■	■	■																						
1-121-86-1	AHN myTip® 200 µL, Gelb		■	■	■	■																						
1-121-C5-1	AHN myTip® 200 µL, Gelb		■	■	■	■																						
1-122-C5-1	AHN myTip® 200 µL, Gelb, steril		■	■	■	■																						
1-127-80-0	AHN myTip® 300 µL, klar		■	■	■	■																						
1-127-C5-0	AHN myTip® 300 µL, klar		■	■	■	■																						
1-128-C5-0	AHN myTip® 300 µL, klar, sterile		■	■	■	■																						
1-130-70-0	AHN myTip® Geloader 200 µL, klar		■	■	■	■																						
1-130-C4-0	AHN myTip® Geloader 200 µL, klar		■	■	■	■																						
1-131-C4-0	AHN myTip® Geloader 200 µL, klar, steril		■	■	■	■																						
1-209-80-2	AHN myTip® 1000 µL, kurz, klar Blau																											
1-202-80-0	AHN myTip® 1000 µL, klar																											
1-202-C4-0	AHN myTip® 1000 µL, klar																											
1-203-C4-0	AHN myTip® 1000 µL, klar, steril																											
1-202-80-2	AHN myTip® 1000 µL, Blau																											
1-202-C4-2	AHN myTip® 1000 µL, Blau																											
1-203-C4-2	AHN myTip® 1000 µL, Blau, steril																											
1-204-80-0	AHN myTip® 1250 µL, klar																											
1-204-C4-0	AHN myTip® 1250 µL, klar																											
1-205-C4-0	AHN myTip® 1250 µL, klar, steril																											
AHN myTip® MaT Makro Pipettenspitzen																												
1-218-61-0	AHN myTip® MaT 5 mL, klar																											
1-222-61-0	AHN myTip® MaT 5 mL, klar, steril																											
1-212-84-0	AHN myTip® MaT 5.5 mL, klar																											
1-216-84-0	AHN myTip® MaT 5.5 mL, klar, steril																											
1-220-50-0	AHN myTip® MaT 10 mL, klar																											
1-224-50-0	AHN myTip® MaT 10 mL, klar, steril																											

■ kompatibel



Kompatibilitätstabelle

Kat.-Nr.	Volumen	CAPP		Sartorius	Brand	Eppendorf	Finnpipette	Gilson Pipetman	Rainin	VWR, Labnet, HTL
		CAPPBravo, ecopipette und CAPPaero Single Channel	CAPP electronic pipettes und CAPPaero Multi Channel	Biohit Eline Proline® Plus mLNE	Transferpette®	Research® plus	Finnpipette Digital	Gilson Pipetman	Rainin Classic, universal shaft (non LTS)	Mechanische Pipetten
AHN myTip® FT Filterspitzen										
2-002-81-0	AHN myTip® FT 10 µL, klar	■	■	■	■	■	■	■	■	■
2-003-C5-0	AHN myTip® FT 10 µL, klar, steril	■	■	■	■	■	■	■	■	■
2-062-81-0	AHN myTip® FT 10 µL XL, klar	■	■	■	■	■	■	■	■	■
2-063-C5-0	AHN myTip® FT 10 µL XL, klar, steril	■	■	■	■	■	■	■	■	■
2-117-80-0	AHN myTip® FT 20 µL, klar	■	■	■	■	■	■	■	■	■
2-118-C5-0	AHN myTip® FT 20 µL, klar, steril	■	■	■	■	■	■	■	■	■
2-119-80-0	AHN myTip® FT 100 µL, klar	■	■	■	■	■	■	■	■	■
2-120-C5-0	AHN myTip® FT 100 µL, klar, steril	■	■	■	■	■	■	■	■	■
2-127-80-0	AHN myTip® FT 200 µL, klar	■	■	■	■	■	■	■	■	■
2-128-C5-0	AHN myTip® FT 200 µL, klar, steril	■	■	■	■	■	■	■	■	■
2-134-80-0	AHN myTip® FT 300 µL, klar	■	■	■	■	■	■	■	■	■
2-135-C5-0	AHN myTip® FT 300 µL, klar, steril	■	■	■	■	■	■	■	■	■
2-202-80-0	AHN myTip® FT 1000 µL, klar	■	■	■	■	■	■	■	■	■
2-203-C4-0	AHN myTip® FT 1000 µL, klar, steril	■	■	■	■	■	■	■	■	■
2-204-80-0	AHN myTip® FT 1250 µL, klar	■	■	■	■	■	■	■	■	■
2-205-C4-0	AHN myTip® FT 1250 µL, klar, steril	■	■	■	■	■	■	■	■	■
AHN myTip® LT Low Retention Pipettenspitzen										
4-002-81-0	AHN myTip® LT 10 µL, klar	■	■	■	■	■	■	■	■	■
4-003-C5-0	AHN myTip® LT 10 µL, klar, steril	■	■	■	■	■	■	■	■	■
4-062-81-0	AHN myTip® LT 10 µL XL, klar	■	■	■	■	■	■	■	■	■
4-063-C5-0	AHN myTip® LT 10 µL XL, klar, steril	■	■	■	■	■	■	■	■	■
4-121-86-0	AHN myTip® LT 200 µL, klar	■	■	■	■	■	■	■	■	■
4-122-C5-0	AHN myTip® LT 200 µL, klar, steril	■	■	■	■	■	■	■	■	■
4-127-80-0	AHN myTip® LT 300 µL, klar	■	■	■	■	■	■	■	■	■
4-128-C5-0	AHN myTip® LT 300 µL, klar, steril	■	■	■	■	■	■	■	■	■
4-202-80-0	AHN myTip® LT 1000 µL, klar	■	■	■	■	■	■	■	■	■
4-203-C4-0	AHN myTip® LT 1000 µL, klar, steril	■	■	■	■	■	■	■	■	■
4-204-80-0	AHN myTip® LT 1250 µL, klar	■	■	■	■	■	■	■	■	■
4-205-C4-0	AHN myTip® LT 1250 µL, klar, steril	■	■	■	■	■	■	■	■	■
AHN myTip® LFT Low Retention Filterspitzen										
5-003-C5-0	AHN myTip® LFT 10 µL, klar, steril	■	■	■	■	■	■	■	■	■
5-063-C5-0	AHN myTip® LFT 10 µL XL, klar, steril	■	■	■	■	■	■	■	■	■
5-118-C5-0	AHN myTip® LFT 20 µL, klar, steril	■	■	■	■	■	■	■	■	■
5-120-C5-0	AHN myTip® LFT 100 µL, klar, steril	■	■	■	■	■	■	■	■	■
5-128-C5-0	AHN myTip® LFT 200 µL, klar, steril	■	■	■	■	■	■	■	■	■
5-135-C5-0	AHN myTip® LFT 300 µL, klar, steril	■	■	■	■	■	■	■	■	■
5-203-C4-0	AHN myTip® LFT 1000 µL, klar, steril	■	■	■	■	■	■	■	■	■
5-205-C4-0	AHN myTip® LFT 1250 µL, klar, steril	■	■	■	■	■	■	■	■	■

■ kompatibel