

Compatibility Chart AHN myPlate® Sealing Membranes with PCR/qPCR 96 Well Plates		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
		Image						
		Standard Well (0.2 mL, height 21 mm)	■	■	■			
		Low Profile Well0 (0.1 mL, height 16 mm)				■	■	■

AHN Sealing Membranes			PCR/qPCR 96 Well Plate						
Cat. No.	Description	Image	AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
3-048-00-0	AHN myPlate® PCR 96 well Sealing Membrane, strong adhesive, operating temperature -40-120°C, PCR use			■	■	■	■	■	■
3-048-01-0	AHN myPlate® qPCR 96 well Optical Sealing Membrane, ultra clear, pressure-sensitive adhesive film, Real-Time PCR (qPCR) use			■	■	■	■	■	■
3-048-02-0	AHN myPlate® qPCR 96 well Optical Sealing Membrane, ultra clear, strong adhesive, operating temperature -40-120°C, Real-Time PCR (qPCR) use			■	■	■	■	■	■
3-048-03-0	AHN myPlate® PCR Aluminium Membrane, strong adhesive, operating temperature -80-120°C, PCR and storage use, piercable			■	■	■	■	■	■