



## Compatibility Chart AHN myPlate® qPCR 96 Well Plate

AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Image						
Standard Well (0.2 mL, height 21 mm)	■	■	■			
Low Profile Well0 (0.1 mL, height 16 mm)				■	■	■

Bio-Rad		qPCR 96 Well Plate						
		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Device	Note							
iCycler IQ			■		■			
iQ5			■		■			
MylQ			■		■			
Chromo4			■		■			
Opticon I & II							■	■
CFX96							■	■

Stratagene/Eppendorf/Biometra		qPCR 96 Well Plate						
		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Device	Note							
Mx4000/3000P/3005P			■		■			
AriaMx							■	■
MasterCycle ep Realplex			■		■			
Optical Thermocycler			■		■			