

Compatibility Chart AHN myTube® qPCR strip 8 tubes

AHN Cat. No.	3-036-03-9	3-036-01-0	3-037-03-0	3-037-04-9
Image				
Standard Well (0.2 mL, height 21 mm)			■	■
Low Profile Well (0.1 mL, height 16 mm)	■	■		

Applied Biosystems (Life Technology)		AHN Cat. No.	3-036-03-9	3-036-01-0	3-037-03-0	3-037-04-9
Device	Holder No.					
7000	4350105				■	■
7300	4367033				■	■
7500	4367033				■	■
7500Fast	4403809		■	■		
StepOne	4375282		■	■		
StepOne Plus	4379983		■	■		
ViiA7 Fast Model (0.1 ml)	4379983		■	■		
ViiA7 (0.2 ml)	4381850				■	■
QuantStudio Fast Model (0.1 ml)	4379983		■	■		
QuantStudio (0.2 ml)	4381850				■	■
QuantStudio 3/5 Fast (0.1 ml)	4379983		■	■		
QuantStudio 3/5 (0.2 ml)	4381850				■	■
3500(xL) (0.1 ml Fast 96 Well)	4410232 (Retainer) / 4410701 (Septa)		■	■		
3500(xL) (0.2 ml Standard 96 Well)	4410230 (Retainer) / 4410701 (Septa)				■	■

Roche		AHN Cat. No.	3-036-03-9	3-036-01-0	3-037-03-0	3-037-04-9
Device	Note					
LightCycler 96			■	■		
LightCycler Nano	Bottom Detection			■		

Compatibility Chart AHN myTube® qPCR strip 8 tubes

AHN Cat. No.	3-036-03-9	3-036-01-0	3-037-03-0	3-037-04-9
Image				
Standard Well (0.2 mL, height 21 mm)			■	■
Low Profile Well (0.1 mL, height 16 mm)	■	■		

Bio-Rad		qPCR strip 8 tubes				
Device	Note	AHN Cat. No.	3-036-03-9	3-036-01-0	3-037-03-0	3-037-04-9
iCycler iQ					■	■
iQ5					■	■
MyiQ					■	■
Chromo4					■	■
Opticon I & II			■	■		
mini Opticon			■	■		
CFX96			■	■		

Stratagene/IT-IS/Open qPCR/Eppendorf/Thermo/Biomtra		qPCR strip 8 tubes				
Device	Note	AHN Cat. No.	3-036-03-9	3-036-01-0	3-037-03-0	3-037-04-9
Mx4000/3000P/3005P					■	■
AriaMx			■	■		
Mygo Pro	Bottom Detection			■		
Open QPCR	Bottom Detection			■		
MasterCycle ep Realplex					■	■
PikoReal 24 Real-Time PCR			■	■		
Toptical Thermocycler					■	■