

# ahn myPlate®


## PCR and qPCR Sealing Membranes





**Usage** - membranes are designed for PCR and qPCR applications and for storage


**Reproducible result** - evaporation rate below 5 %

**Ultra clear membranes** - the optical sealing membranes 3-048-01-0 and 3-048-02-0 show consistent and high light-transmitting rates for stable data



## Ordering information

Volume	Features	Pack Type	Sales Unit	Cat. No.
<b>96-well Plate Sealing Membrane</b> 	<ul style="list-style-type: none"> <li>operating temperature -40 ~120°C</li> <li>PCR use</li> <li>strong adhesive</li> </ul>	bag	10x10	3-048-00-0
<b>96-well Plate Optical Sealing Membrane</b> 	<ul style="list-style-type: none"> <li>operating temperature -40 ~120°C</li> <li>Real-Time PCR(qPCR) use</li> <li>strong pressure sensitive adhesive</li> </ul>	bag	10x10	3-048-01-0
<b>96-well Plate Optical Sealing Membrane</b> 	<ul style="list-style-type: none"> <li>operating temperature -40 ~120°C</li> <li>Real-Time PCR(qPCR) use</li> <li>strong adhesive</li> </ul>	bag	10x10	3-048-02-0
<b>Aluminium Membrane</b> 	<ul style="list-style-type: none"> <li>operating temperature -80 ~120°C</li> <li>PCR or Storage use</li> <li>strong adhesive</li> <li>piercable</li> </ul>	bag	10x10	3-048-03-0

Accessory	Features	Pack Type	Sales Unit	Cat. No.
<b>Decapping/ Sealing Tool for PCR Strips PCR Membranes</b> 	<ul style="list-style-type: none"> <li>PCR strip decapping: 3-037-01-0, 3-037-02-0, 3-037-03-0, 3-037-04-9</li> <li>Membrane sealing: 3-048-00-0, 3-048-01-0, 3-048-02-0, 3-048-03-0</li> </ul>	bag	1	3-000-BOP-0

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com

www.ahn-bio.de

