

ahn myLab CLC-02

Clinical Centrifuge 4500 rpm

High-tech microprocessor - Advanced microprocessor technology ensures dependable accuracy and precision with each spin cycle

Intuitive user controls - User-facing control panel provides quick access to all device features for quick and easy programming of the clinical centrifuge

Maintenance-free DC motor - Brushless DC motor with variable speeds provides reliable performance with minimum maintenance or downtime

Combi rotor - Convenient 8-slot rotor design that fits 8 x 15 mL or 4 x 50 mL tubes at maximum capacity

Excellent safety features - Low noise, imbalance detection, and lid-lock safety mechanisms help to maintain an accident-free laboratory environment

Compact and portable - Space-saving centrifuge design takes up a small amount of workbench real estate while a lightweight chassis makes it easy to move from one point to another



Ordering information

Specifications	
Display	Large, easy-to-read digital display
Motor type	Brushless DC motor
Rotor	Combi rotor
Rotor capacity	8 x 15 mL / 4 x 50 mL tubes
Relative Centrifugal Force	1985 x g (max.)
Speed settings	500 – 4500 rpm (max.)
Timer setting	1 – 99 mins and infinity
Noise level	< 60 dB
Device weight	4 Kgs (approx.)
Dimensions (WxDxH)	330 x 265 x 161 mm
Lid-lock	Yes
Imbalance detection	Yes

Description	Cat. No.
AHN myLab® CLC-02 Clinical Centrifuge 4500 rpm	7-011-00-0

Delivery Package
1 pc. Clinical centrifuge 4500 rpm
1 pc. Power supply adaptor
1 pc. T-Allen Key
1 pc. 8-slots rotor (pre-installed)
4 pcs. 15mL Reduction Adaptors
4 pcs. 15mL Caps
8 pcs. Reduction Adaptors & Rubber Caps
1 pc. Instruction manual
1 pc. Warranty card

ahn



+49(0)3631/65242-0



+49(0)3631/65242-90



Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com



www.ahn-bio.de

