

# ahn myPlate®

## Petri Dishes

**Superior quality, contaminant-free petri dishes** - manufactured in a cleanroom facility, resulting in sterile dishes that are free of endotoxins, pyrogens and other common contaminants.

**Straightforward packaging** - supplied in simple packet options. Each pack has 10 or 20 dishes per bag, allowing users to choose the most optimal quantity for their media preparation. This reduces the risk of wasted media and contaminated dishes.

**Space-saving** - with a standard sized diameter of 90mm, AHN Petri Dishes have the perfect balance between plating area and space-saving dimensions. Additionally, these dishes benefit from being both lightweight and fully stackable.

**Excellent airflow** - AHN Petri Dishes are supplied in vented and non-vented formats to accommodate different microbial culture needs. Vented plates have a lip that allows for airflow and moisture control whilst non-vented dishes limit airflow and maximise moisture retention.



## Ordering information

Description	Diameter	Height	Vented	Sterile	Aseptic	Colour	Pack Type	Sales Unit	Cat. No.
AHN myPlate® Petri Dish	91.5mm	14.35mm			ISO class 6	clear	bag	48x10	3-510-48-0
AHN myPlate® Petri Dish	91.5mm	14.35mm		Yes	ISO class 6	clear	bag	48x10	3-511-48-0
AHN myPlate® Petri Dish	91.5mm	14.55mm	Yes, 3 cams		ISO class 6	clear	bag	48x10	3-512-48-0
AHN myPlate® Petri Dish	91.5mm	14.55mm	Yes, 3 cams	Yes	ISO class 6	clear	bag	48x10	3-513-48-0
AHN myPlate® Petri Dish	91.5mm	14.35mm			ISO class 6	clear	bag	24x20	3-520-48-0
AHN myPlate® Petri Dish	91.5mm	14.35mm		Yes	ISO class 6	clear	bag	24x20	3-521-48-0
AHN myPlate® Petri Dish	91.5mm	14.55mm	Yes, 3 cams		ISO class 6	clear	bag	24x20	3-522-48-0
AHN myPlate® Petri Dish	91.5mm	14.55mm	Yes, 3 cams	Yes	ISO class 6	clear	bag	24x20	3-523-48-0

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com

www.ahn-bio.de

