## ahn reagent Reservoirs

## 10 mL | 25 mL | 50 mL

Construction - made out of a unique blend of Polystyrene

Extra large, sturdy base with broad edges - enable it to stand firmly on your workbench

Graduation marks - for easy, accurate dispensing / helps identify any remaining liquid

Angled interior supports the collection of the liquid and decreases reagent waste

Spillproof - no need to worry about spillage while you are processing liquids

Easy pouring - to store any remaining liquid you can empty the reservoir using the corners

Compact storage - the reservoirs are easy to stack and the lip at the bottom of the reservoir allows them to be easily separated from each other



## Ideal for using with all AHN Single and Multi channel Pipettes

- to suit the various pipettes our reservoirs come in three different sizes:
- 10 mL, V-shaped bottom (ideal with 8- and 16-channel pipettes)
- 25 mL, flat bottom (ideal with multi channel pipettes)
- 50 mL, V-shaped bottom (ideal with 12-channel pipettes)



## **Ordering information**

Description	Pack Type	Sales Unit	Cat. No.
AHN Reagent Reservoirs 10 mL	bag	5x100	3-701-11-9
AHN Reagent Reservoirs 25 mL	bag	5x100	3-702-11-9
AHN Reagent Reservoirs 50 mL	bag	5x100	3-703-11-9











